



**icon**  
**are**  
smart independence



# Product Catalogue

EDITION 3.0

# Ta-Da® Chair

## A Seat Stick

Walking Stick Converts to Chair Automatically



Water Proof



Comfortable



Safety Lock



Double Locking



Auto Leg-unlocking



Auto Leg-locking

The Ta Da® chair from Iconicare is a brilliant lightweight walking stick that converts into a tripod seat at the touch of a button.

Made from strong lightweight aluminium with a solid plastic seat, the Ta Da® chair only weighs 920 grams. The three feet with rubber Ferrules help keep your customer steady when walking. When your customer needs a break, they can stop where they are.

It is ideal for use in the countryside, walking around the garden and in town. Comes with a high quality luxury carry bag with a strap to put on over the shoulder.

The Iconicare TaDa® Chair is of the highest quality and comes with a money back guarantee if your customers are not completely satisfied.

### Features

- Ta-Da®, is a versatile mobility aid that simultaneously provides the user with a walking stick for support and a seat for resting.
- This design, which instantly converts a cane into a chair, offers rest for tired legs and knees at a moment's notice.
- With its unique design, the seat can adjust itself according to the angle of the user's body to maximize sitting comfort.
- Once converted into a stick, the bulk weight of the collapsed portable chair is up near the stick's handle, thereby reducing the user's effort to swing and carry the stick.
- Several ground breaking patented designs make the product safe and easy to use. Users will be able to experience its extraordinary functionality and practicality.
- Seat height 51cm (20")
- Weight Capacity: 100kgs.
- Comes with a high quality luxury carry bag with a strap to put on your shoulder.
- Designed for people who are aged 16 or above are able to walk independently and have good sitting balance while sitting.

# Shoe Helper

Shoe Horn & Shoe Remover all in one.



Our Shoe Helper is 2.5 feet (30" / 76.5cm) tall. One end of the Shoe Helper is a traditional shoe horn to help get your shoes on. The other end is a horseshoe shaped aid to lock shoes in place so the shoe can be taken off easily.

Our users include the elderly, people with back problems and limited mobility and pregnant ladies

Iconicare have an exclusive distribution arrangement for the UK.

## Features

- Make putting shoes on and taking off an easy task
- Take shoes on or off simply and quickly
- The Shoe Helper is lightweight, durable and free standing
- Won't slip or damage footwear
- Available in wood (two colours) and plastic
- Simple and stylish.
- Like all our products, the Iconicare Shoe Helper is of the highest quality and comes with a money back guarantee if your customers are not completely satisfied.



Iconicare specialises in unique, individual retail products giving smart independence

We offer your customers the ability to purchase these products supplied directly from our UK warehouse at Trade prices.

For further information on these products please request our separate Iconicare brochure at [www.iconicare.co.uk](http://www.iconicare.co.uk) or contact [sales@iconicare.co.uk](mailto:sales@iconicare.co.uk) for further information.

## Terms & Conditions

Title of goods changes hands at customer premises.

Urathon UK is the Consignee/Importer.

Urathon UK is responsible for organising and payment of Sea freight, Marine Transit Insurance, all duties, taxes, handling, clearance and haulage fees.

PAYMENT TERMS: Payment by debit or Credit card with order

DELIVERY: £500.00 Carriage paid, Carriage for Individual items chargeable at £5.50 per item.

## Iconicare Smart Independence

Customer Services:

**+44 (0) 7366568767**

**[sales@iconicare.co.uk](mailto:sales@iconicare.co.uk)**

Thane House  
Hilmarton  
Calne  
SN11 8SB



**[www.iconicare.co.uk](http://www.iconicare.co.uk)**