

# COMMUNITY EQUIPMENT CATALOGUE

Mobility & Rehabilitation



Call | 01249 760581  
Visit | [www.urathon.com](http://www.urathon.com)  
Email | [sales@urathon.com](mailto:sales@urathon.com)



# **Urathon**<sup>®</sup>

Specialist suppliers to the mobility & rehabilitation industry

# Urathon®

For nearly 30 years, Urathon has been a trusted name in the healthcare and mobility sector, supplying high-quality products to many of the UK's leading providers. We support some of the country's largest providers, including NHS services, hospitals, care homes, and mobility retailers, by delivering reliable, high-quality solutions.

We offer a complete supply chain solution, from product design and development to manufacturing, branding, and distribution. Our operations are managed from our headquarters in Wiltshire, UK, where our warehouse holds a wide range of ready-to-ship healthcare and mobility items. For larger or custom orders, we also provide direct container shipments from our wholly owned, ISO-certified manufacturing facility in China.

Urathon specialises in Aids for Daily Living (ADLs), mobility accessories, and medical equipment, with a strong focus on quality, compliance, and user comfort. Whether you need off-the-shelf products or bespoke branded solutions, our UK-based team ensures smooth coordination, reliable delivery, and responsive customer support.

With Urathon, you're not just buying products, you're partnering with a team committed to helping you deliver better care, streamline your operations, and meet the evolving needs of the healthcare sector with confidence and efficiency.



## UK STOCKED


Distributed from our UK warehouse



## FROM CHINA BY MOQ

Available via direct mixed container, or to our UK facility. Priced upon request.

# CONTENTS





### 4 Wheeled Aluminium Rollator

- Foldable
- 200mm PU Castor
- Seat, Back Rest and Bag
- Length: 670mm
- Width: 630mm
- Width between arms: 445mm
- Height: 800 - 980mm
- Seat height: 550mm
- Max weight: 165kg
- BS EN ISO 10993-10



90387



### 3 Wheeled Aluminium Rollator

- Lightweight Grey powder coated Aluminium frame
- Easy to fold
- Height adjustable push handles
- Push to lock brakes
- 200mm PU castor wheels
- Bag included
- Height: 800 - 980mm
- Length: 670mm
- Total Width at rear: 685mm
- Max weight: 125kg
- BS EN ISO 10993-10



90386



## Moving & Handling

Discover our 'Moving and Handling' range, a dedicated collection of products designed to assist occupational therapists and caregivers in providing safe and efficient patient care.

Within this section, you'll find a comprehensive range of equipment and aids specially crafted to facilitate the seamless movement of patients from one place to another. Our selection includes items that prioritize patient comfort, safety, and ease of use, enabling professionals to deliver top-quality care.

With these valuable tools at your disposal, you can confidently and effectively manage patient transfers and handling, ensuring their well-being and comfort at every step.

### RS1 Patient Turner

- Enables sit-to stand for semi-ambulant users
- Height adjustable handle
- Height adjustable knee pads
- Slip resistant turning disc
- Easy to use and manoeuvre
- Length: 469mm
- Width: 445mm
- Height: 786 - 1215mm
- Castor Size: 75mm
- Max weight: 200kg
- BS EN ISO 10993-10



90330



### RS2 Patient Transfer Assist

- Low step on height base for ease with clients who have minimal leg movement
- Anti slip base ensures the client is safe and secure in transfer
- Low friction castor for ease of manoeuvrability
- One foot touch brake keeps the base stable when embarking and disembarking from the raiser
- Versatile handle design ensures carers of all shapes and sizes can easily manoeuvre it
- Adjustable knee pad gives the client support and comfort when in use
- L: 635mm, W: 607mm, H: 1186mm
- Max weight: 160kg
- BS EN ISO 10993-10



90329



### Ladder Belt (Sliding Support Strap)

- For use on RS2 Patient Transfer Assist
- Length: 2000mm
- Width: 100mm
- Height: 20mm
- Max weight: 160kg



90359



### Patient Handling Belt

- For use on RS2 Patient Transfer
- 3 sizes available: Small, Medium and Large
- Length (mm): S:1050 / M:1300 / L:1500
- Width (mm): S:975 / M:1150 / L:1350
- Height: 100mm
- Max weight: 160kg



Small	90333S
Medium	90333M
Large	90333L

## RS4 Patient Transfer Assist

- Light Weight, Strong and Robust frame
- Easy flip down support pads
- Easy to open leg mechanism. Works in conjunction with most types of seating and chairs
- Enables users to participate in transfers, maintaining patient independence
- Ideal for those users who can raise from a seated position to a standing position
- Adjustable kneepad
- Compact narrow frame. Manoeuvre in small spaces
- Length: 820mm, Width: 530-730mm, Height: 1060mm
- Seat height: 620-670mm
- Max weight: 185kg
- BS EN ISO 10993-10



3" Subassembly  
90580



3" Full Assembly  
90720



4" Subassembly  
90724



4" Full Assembly  
90733

# AQUA LIBRA

## STAINLESS STEEL SHOWER CHAIR RANGE



CONTACT OUR SALES TEAM ON  
**01249 760581**  
sales@urathon.com

# Aqua Libra AP

ATTENDANT PROPELLED

# Aqua Libra SP

SELF PROPELLED



3 HANDLE OPTIONS AVAILABLE



	Straight Handle		Push Handle		Pram Handle	
	AP	SP	AP	SP	AP	SP
17"	90646-50001	90649-50001	90646-50010	90649-50010	90646-50023	90649-50023
19"	90647-50002	90650-50002	90647-50011	90650-50011	90647-50024	90650-50024
21"	90648-50003	90651-50003	90648-50012	90651-50012	90648-50025	90651-50025



## FEATURES AND BENEFITS

- Fully modular construction – ability to adapt the chair to suit any user
- Stainless steel, white powder coated frame – cannot rust and lends itself well to constant recycling if required
- Height adjustment frame as standard, designed to be able to fit over most standard toilets
- Large commode pan capacity
- 17", 19" and 21" frame width options
- Available in both attendant and self propelled variants
- Wide range of seat variances – ideal for all toileting requirements
- Adjustable / easily removable arms can be supplied with various width options
- Easily removable backrest with 3 handle options available
- Easy to clean seat and backrest materials
- Fully removable swing away footrests with simple length adjustment and nylon footplates
- Max weight: 160kg



### AP Aluminium Shower Commode Chair

- Removable PU Back and PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 460mm
- Seat height: 510 - 560mm
- Height: 990 - 1010mm
- Length: 900mm
- Width: 570mm
- BS EN ISO 10993-10



90209



### Self Propul Wheel Kit

- Turn your AP into a SP Aluminium Shower Commode Chair
- Removable PU Back and Horseshoe Seat with Aperture Plug
- 2 x Braked Nylon Castors and 2 x 24" MAG Wheels with PU Tyres and Q.R Pins
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 460mm
- Seat height: 510 - 560mm
- Height: 990 - 1010mm
- Length: 1085mm
- Width: 660mm
- BS EN ISO 10993-10



90351



### AP HD Aluminium Shower Commode Chair

- Removable PU Back  
PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 660mm
- Seat height: 510 - 560mm
- Height: 990 - 1010mm
- Length: 900mm
- Width: 800mm
- BS EN ISO 10993-10



90350



## AP XXL Aluminium Shower Commode Chair

- Removable PU Back, PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 700mm
- Seat height: 510 - 535mm
- Length: 900mm
- Width: 850mm
- Height: 985-1010mm
- Max weight: 300kg
- BS EN ISO 10993-10



90513



## AP Aluminium Shower Commode Chair, Tilt in Space

- Lightweight aluminium frame - Easy to manoeuvre and offers long-term rust-free durability.
- Tilt and space chair – ideal for showering people comfortably
- 4 x Large, braked, robust castors for durability and manoeuvrability.
- Polyurethane seat, seat-infill, back and armrests are water-resistant, easy to clean and comfortable.
- Full height adjustable to suit all heights of the users
- Simple swing away and height adjustable footrests
- Easy to clean seat and removable back
- Flip up arm rests – For easy side transfer
- Large capacity commode pan and lid
- BS EN ISO 10993-10



AVAILABLE NOW  
Tension adjustable mesh backrests

90581



**Height Adjustable Chrome Plated Glideabout Commode**

- Height Adjustable
- Height: 985 - 1085mm
- Seat height: 525 - 625mm
- Max weight: 160kg



90294



**Chrome Plated Glideabout Commode**

- 4 x Castors
- Removable, Grey Aperture Seat with Padded Vanity Cover
- Removable Arm Rests, Backrest, Leg Rests
- Potty
- Grey Vinyl Upholstery
- Length: 840mm
- Width: 525mm
- Height: 985mm
- Seat height: 535mm
- Max weight: 190kg



90099



**Mild Steel Glideabout Commode**

- 4 x Castors
- Removable, Grey Aperture Seat with Padded Vanity Cover
- Removable Arm Rests, Backrest, Leg Rests
- Potty
- Grey Vinyl Upholstery
- Length: 922mm
- Width: 530mm
- Height: 962mm
- Seat height: 528mm
- Max weight: 190kg
- BS EN ISO 10993-10



90517



**Height Adjustable Mild Steel Glideabout Commode**

- Height adjustable
- Height: 985 - 1085mm
- Max weight: 165kg
- BS EN ISO 10993-10



90675



**KD Glideabout Commode**

- Easy to Disassemble
- Removable Aperture
- Padded Vinyl Seat
- Removable Back, Arms, Swing Angle Leg Rests
- Length: 930mm
- Width: 540mm
- Height: 980mm
- Seat height: 500mm
- Max weight: 165kg
- BS EN ISO 10993-10



90519



### Shower Stool

- Height Adjustable
- Powder Coated Steel Frame
- Easy to remove PP Seat
- Aluminium option
- Length: 370 - 395mm
- Width: 495 - 540mm
- Height: 560 - 710mm
- Seat height: 480 - 630mm
- Max weight: 160kg
- BS EN ISO 10993-10



90115



### Drop Down Rail with Support Leg

- Length: 692mm
- Width: 100mm
- Height: 750mm
- Max weight: 190kg
- BS EN ISO 10993-10



90146



2"

### Raised Toilet Seat, No Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90635



4"

### Raised Toilet Seat, No Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90636



6"

### Raised Toilet Seat, No Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90637



2"

### Raised Toilet Seat, With Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90660



4"

### Raised Toilet Seat, With Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90661



6"

### Raised Toilet Seat, With Lid

- Ergonomic design, strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg



90662



**Toilet Surround Height and Width Adjustable**

- Height & Width Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 545 - 670mm
- Height: 645 - 795mm
- Max weight: 190kg
- BS EN ISO 10993-10



90438



**Toilet Surround**

- Height Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 524mm
- Height: 645 - 795mm
- Max weight: 190kg
- BS EN ISO 10993-10



90028



**Static Shower Chair**

- Height Adjustable
- Length: 495mm
- Width: 620 - 645mm
- Height: 775 - 925mm
- Seat height: 420 - 570mm
- Max weight: 165kg
- BS EN ISO 10993-10



90300



**4 Part Bath Step**

- Length: 450mm
- Width: 350mm
- Height: 200mm
- Max weight: 500kg



90116



**Male Urinal Bottle with lid**

- Length: 126mm
- Width: 98mm
- Height: 225mm



90597



**Bath Seat**

- Strong suction feet
- Available in 150mm, 250mm and 350mm height options
- Length: 280mm
- Width: 460mm
- Max weight: 190kg



90040 6"  
90041 8"  
90042 12"



### Toilet Frame, Height Adjustable

- Injection Moulded Seat
- L: 443 - 471mm
- W: 607 - 612mm
- H: 622 - 771mm
- Seat H: 428 - 578mm
- Max weight: 190kg
- BS EN ISO 10993-10



90341



### Toilet Frame, Width & Height Adjustable

- Injection Moulded Seat
- L: 440 - 468mm
- W: 570 - 720mm
- H: 622 - 771mm
- Seat H: 428 - 578mm
- Max weight: 190kg
- BS EN ISO 10993-10



90532



### Swivelling Bath Chair

- Adjustable Width
- Width: 770 - 820mm
- Height: 530mm
- Width between arms: 600 - 700mm
- Max weight: 160kg
- BS EN ISO 10993-10



90071



### Toilet Frame, Clip on Seat & Chute

- Height Adjustable
- Length: 465 - 495mm
- Width: 558mm
- Height: 690 - 835mm
- Seat height: 483 - 628mm
- Max weight: 190kg
- BS EN ISO 10993-10

Version with lid available (90698)



90120



### 26" Bath Board 4 Slat

- Adjustable securing brackets
- Depth: 225mm
- Length: 660mm
- Height: 110mm



90210



### 27" Bath Board 4 Slat

- Adjustable securing brackets
- Depth: 225mm
- Length: 685mm
- Height: 110mm



90211



### 28" Bath Board 4 Slat

- Adjustable securing brackets
- Depth: 225mm
- Length: 710mm
- Height: 110mm



90212



### 26" Bath Board 6 Slat

- Adjustable securing brackets
- Depth: 360mm
- Length: 660mm
- Height: 200mm
- Max weight: 190kg



90254



### 27" Bath Board 6 Slat

- Adjustable securing brackets
- Depth: 360mm
- Length: 690mm
- Height: 200mm
- Max weight: 190kg



90255



### 28" Bath Board 6 Slat

- Adjustable securing brackets
- Depth: 360mm
- Length: 710mm
- Height: 200mm
- Max weight: 190kg



90256

## Bedroom Aids

Welcome to our 'Bedroom Aids' section, where you'll find a curated selection of essential products designed to enhance safety, comfort, and convenience in your bedroom. Explore our range of bed rails, bed supports, and commodes, thoughtfully chosen to provide the support you need to create a secure and comfortable bedroom environment. Whether you're looking for assistance with mobility or seeking solutions to improve your overall quality of life, our bedroom aids are here to help you rest easy.



### Double Bedstick

- Adjustable Width
- Chrome Plated
- Length: 335mm
- Width: 960 - 1555mm
- Height: 440mm
- Max weight: 140kg



90355



### Double Bedstick

- Adjustable Width
- Powder Coated White
- L: 335mm W: 960 - 1555mm H: 440mm
- Max weight: 140kg
- BS EN ISO 10993-10



90354



### Bedstick Slatted Bed

- L: 830mm W: 625mm H: 470mm
- Max weight: 190kg
- BS EN ISO 10993-10



90321



### Bedstick For Divan Bed

- Length: 810mm
- Width: 625mm
- Height: 470mm
- Max weight: 190kg
- BS EN ISO 10993-10



90322



### Straps (PAIR) for UEL90321/90322

- Length: 2500mm
- Width: 25mm
- Max weight: 50kg



90324



### Bed Grab Rail

- Length: 620mm
- Width: 365mm
- Height: 470mm
- Max weight: 160kg
- BS EN ISO 10993-10

90521



### Bed Lever with Straps

- Length: 620mm
- Width: 365mm
- Height: 470mm
- Max weight: 160kg
- BS EN ISO 10993-10

90325



### Stacking Commode Chair

- Height Adjustable
- Fixed Arms
- HF Vinyl Welded Padded Seat & Back
- Brown Powder Coated Steel Frame
- Anti Rattle Ring & DIP Moulded Ferrule
- Commode Pan & Lid
- Length: 520 - 545mm
- Width: 551mm
- Height: 785 - 930mm
- Seat Height: 470 - 620mm
- Max weight: 180kg
- BS EN ISO 10993-10



90096



### Stacking Commode Chair with Removable Arms

- Height Adjustable
- Removable Arms
- HF Welded Padded Seat & Back
- Brown Padding and Brown Powder Coated Steel
- Anti Rattle Ring & DIP Moulded Ferrule
- Commode Pan & Lid
- Length: 460 - 505mm
- Width: 550mm
- Height: 780 - 930mm
- Seat Height: 490 - 640mm
- Max weight: 180kg
- BS EN ISO 10993-10



90323



### Bed Bumpers with Mesh Divide

- Blue with black mesh
- Length: 2000mm
- Width: 710mm
- Height: 10mm



90522



### Bed Bumpers with Mesh Divide

- Cream with white mesh
- Length: 2000mm
- Width: 710mm
- Height: 10mm



90560



### Postural Wedge

- Non slip material
- Length: 530mm
- Width: 275mm
- Height: 210mm



90397

## Hospital Portering Chair

- Full PU Seat and back – Durable, robust, easy to clean and reduces any cross-infection risk. Ideal in heavy usage situations
- Strong balanced footplate – Durable with the ability to take the full user weight without tipping
- The footplate is made from HDPE plastic – Robust and ideal for infection control. Competitors tend to use wood which is a cross infection risk
- Stainless steel frame – Ensure product longevity. Allows the product to be used outside without fear of rust. Easy to wipe surface reduces risk of cross infection. Competitors are made from mild steel that can rust when exposed to the inclement weather
- Max user weight of 190kg
- Moulded drop-down arms – Allows for quick and safe side transfers
- Accessory mounts – Allows IV poles etc to be used with the chair
- Xray/notes pocket as standard – Enables the users records to be transported safely with the user from department to department.
- BS EN ISO 10993-10
- Size: 600 x 865 x 1040mm
- Seat Width: 540mm
- Seat Depth: 410mm
- Seat Height: 540mm
- Backrest Height: 480mm



90337

## Around The Home

In our 'Around the Home' range you'll find a thoughtfully curated collection of products designed to improve the functionality and comfort of your living space. Explore our range of footstools, overbed tables, and perching stools, each carefully selected to enhance your daily life. Whether you're seeking extra support while sitting or assistance with daily tasks, our home-focused aids are here to make your surroundings more accommodating and user-friendly. Elevate your home environment and enjoy increased convenience and ease with our selection of purposeful products.



### Overbed Table

- Height Adjustable
- NO Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Top
- Length: 697mm
- Width: 413mm
- Height: 630 - 930mm
- Max weight: 10kg
- BS EN ISO 10993-10



90290



### Overchair Table

This curved overchair/bed table is both width and height adjustable allowing the table to be set up perfectly to the requirements of the user. This table is great for use with chairs, single beds and wheelchairs.

The kidney shaped top is designed to curve gently around the user reducing the need for them to reach forward, it also allows them to rest their arms comfortably on the table top whilst still having a large usable work area.

- Height Adjustable
- Width: 1010 - 1160mm
- Depth: 510mm
- Height: 710-910mm
- Max Weight: 15kg



90615



### Overbed Table with Castors

- Height Adjustable
- 4 x 2" Braked Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Table Top
- Length: 680mm
- Width: 410mm
- Height: 680 - 980mm
- BS EN ISO 10993-10



90291



### Vinyl Footstool

- Height Adjustable
- 4 x Ferrule
- Length: 500mm
- Width: 280mm
- Height: 370 - 520mm
- Seat height: 380 - 535mm
- Max weight: 100kg
- BS EN ISO 10993-10



90454



### Vinyl Footstool on Castors

- Height Adjustable
- 4 x 40mm Castors
- Length: 500mm
- Width: 280mm
- Height: 370 - 520mm
- Seat height: 375 - 525mm
- Max weight: 100kg
- BS EN ISO 10993-10



90455



### Perching Stool

- HF Vinyl Padded Seat
- Adjustable Height
- Length: 415- 460mm
- Width: 561 - 620mm
- Height: 540 - 686mm
- Seat height: 540 - 686mm
- Max weight: 165kg
- BS EN ISO 10993-10



90167



### Perching Stool with Arms

- HF Padded Seat
- Steel Loop Arms, No Back
- Adjustable Height
- Length: 415- 460mm
- Width: 561 - 620mm
- Height: 723 - 868mm
- Seat height: 540 - 686mm
- Max weight: 165kg
- BS EN ISO 10993-10



90168



### Perching Stool with Arms and Back

- HF Vinyl Seat
- Steel Loop Arms & Steel Back
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 825 - 965mm
- Seat height: 540 - 685mm
- Max weight: 165kg
- BS EN ISO 10993-10



90169



### Perching Stool with Arms and HF Back

- HF Vinyl Seat
- Steel Loop Arms & HF Back
- Adjustable Height
- Length: 440 - 460mm
- Width: 560 - 620mm
- Height: 830 - 975mm
- Seat height: 540 - 685mm
- Max weight: 165kg
- BS EN ISO 10993-10



90245



### Bariatric Perching Stool with Arms & HF Back

- HF Seat and HF Back
- Adjustable Height
- Length: 500 - 550mm
- Width: 690 - 750mm
- Height: 950 - 1095mm
- Seat height: 570 - 715mm
- Max weight: 254kg
- BS EN ISO 10993-10



90327



### Perching Stool with PU Seat

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 535 - 685mm
- Seat height: 535 - 680mm
- Max weight: 165kg
- BS EN ISO 10993-10



90270



### Perching Stool with Large Blue PU Seat and Steel Arms

- Adjustable Height
- Length: 420 - 465mm
- Width: 560 - 635mm
- Height: 722 - 862mm
- Seat height: 535 - 685mm
- Max weight: 165kg
- BS EN ISO 10993-10



90530



### Perching Stool with PU Seat, Steel Arms and Back

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 635mm
- Height: 825 - 965mm
- Seat height: 535 - 685mm
- Seat width: 368mm
- Seat depth: 278mm
- Max weight: 165kg
- BS EN ISO 10993-10



90272



### Perching Stool with PU Seat, Steel Arms and Padded Back

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 635mm
- Height: 865 - 1005mm
- Seat height: 535 - 685mm
- Seat width: 368mm
- Seat depth: 278mm
- Max weight: 165kg
- BS EN ISO 10993-10



90273



### Bariatric Perching Stool with PU Seat, PU Back & Arms

- Adjustable Height
- Length: 500 - 550mm
- Width: 690 - 750mm
- Height: 950 - 1095mm
- Seat height: 570 - 745mm
- Max weight: 280kg
- BS EN ISO 10993-10



90346



### KD Perching Stool

- Removable Arms, Legs and Back
- Steel Arms
- Height Adjustable
- Max weight: 160kg
- BS EN ISO 10993-10



90626



### Two Tier Trolley with Brakes

- Push to lock brakes
- Length: 530mm
- Width: 585mm
- Height: 762 - 912mm
- Castor size: 150mm
- Max weight: 125kg
- BS EN ISO 10993-10



90515



### 2 Tier KD Dinner Trolley

- Height Adjustable
- 2 x Removable Shelves
- Length: 501mm
- Width: 445mm
- Height: 845 - 945mm
- Distance between trays: 266mm
- Max weight: 180kg
- BS EN ISO 10993-10



90072



### Half Step

- Adjustable Feet
- Removable slip resistant mat
- Plastic - will not rust
- Length: 430mm
- Width: 735mm
- Height: 112mm
- Max weight: 190kg



90141



### Back Rest

- Folds flat for easy storage
- Sturdy frame with breathable mesh fabric
- Adjustable soft pillow
- 5 levels of angle adjustment
- Width: 580mm
- Depth: 515mm
- Height: 370 - 650mm
- Max weight: 127kg



90668



# ***Urathon***<sup>®</sup>

## **UK WAREHOUSING**

Our dedicated UK warehouse allows for quick dispatch to you. No overseas shipping delays. Receive your products quickly, on time and when you need them. That's our promise to you.

## **SUPPLY FROM STOCK**

We fulfil product orders quickly from stock within our self managed UK warehouse. This allows us to always provide fast and reliable turnaround times with next day delivery available.

## **OUTSTANDING CUSTOMER CARE**

Customer care is at the heart of our business. Your satisfaction is paramount and we always aim to go above and beyond.

Call **01249 760581**

Email **sales@urathon.com**

Visit **www.urathon.com**

***Urathon***<sup>®</sup>

**SPECIALIST SUPPLIERS TO THE MOBILITY  
& REHABILITATION INDUSTRY**



Urathon Ltd. Company No. 02766713