

Telephone **01249 760581** or email **sales@urathon.com**

2026

Product Catalogue




Urathon[®]

SPECIALIST SUPPLIERS TO THE MOBILITY
& REHABILITATION INDUSTRY

Urathon®

Urathon have been designing and brand manufacturing for many of the healthcare industries leading suppliers for almost 30 years including many of the Uks largest providers to Nursing homes, Hospitals, NHS community equipment stores and Mobility retailers.

We offer a complete seamless supply chain solution for all your Healthcare and Mobility product requirements either from our large fully stocked warehousing in Wiltshire, UK or by full mixed containers under your own brand directly from our wholly owned manufacturing base in China.



CONTENTS

THOR

ULTRA - LIGHT FOLDING SCOOTER

Discover effortless mobility with this ultra-light folding scooter, combining sleek modern design with a full carbon-fibre easy-fold frame weighing just 13kg. Exceptionally portable and simple to store, it's ideal for everyday use and travel.

Comfort features include padded upholstery and fold-up, removable armrests, while clear, easy-to-use controls offer horn, speed adjustment, and indicators for confident riding. Quiet 180W motors provide smooth performance and a range of up to 15.5 miles.

A lightweight, removable aircraft-safe lithium battery, USB-C charging, and an optional magnetic phone holder complete this smart, convenient mobility solution—designed to keep you moving with ease.



Lightweight Carbon Fibre Finish

90760

HERO

CARBON FIBRE POWERCHAIR

Revolutionary frame design manufactured from aviation grade lightweight carbon fibre, available in stunning Ruby red and Matt Onyx black.

Easy to fold and store, luxury padded backrest and memory foam cushion.

Advanced, intuitive easy to use controls, fold up armrests for easy access.

2x 250w powerful brushless motors powered by an easy removable aircraft-safe lithium battery gives an impressive 10 mile range, capable of climbing 12 degree slopes with 7cm ground clearance.



Matt Carbon Fibre Finish





LEO Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Breathable fabric seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large 10" Front and 8" Rear wheels with PU Flat Free Tyres
- High quality easily removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Max Weight: 130kg

| | |
|--------|-------|
| Small | 90653 |
| Medium | 90652 |



ARIES Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Uprated soft padded seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large uprated 10" Front and 8" Rear wheels with PU Flat Free Tyres
- High quality easily removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Seat Height: Small 540mm, Medium 620mm
- Max Weight: 130kg

| | |
|--------|-------|
| Small | 90655 |
| Medium | 90654 |



TAURUS Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Luxury Genuine stitched leather finishing on padded seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large uprated 10" Front and 8" Rear wheels with chunky treaded PU Flat Free Tyres
- Genuine stitched leather removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Seat Height: Small 540mm, Medium 620mm
- Max Weight: 130kg



Stitched Leather

| | |
|--------|-------|
| Small | 90657 |
| Medium | 90656 |

CIRRUS

LIGHTWEIGHT 4 WHEEL ROLLATOR

Discover the perfect blend of elegance and practicality with our Cirrus Collection of foldable aluminium rollators.

Lightweight Aluminium Frame: Engineered from high-quality aluminium, our rollators offer exceptional strength and durability while remaining lightweight for easy manoeuvrability.

Smooth Rolling Wheels: Glide effortlessly on our PP+PVC 8" wheels, designed to provide a smooth and stable ride on any surface.

Sleek Metallic Finish: Make a statement with our rollators' sleek metallic finish, combining sophistication with practicality.

Compact Design: With dimensions tailored for easy navigation, our rollators offer the perfect balance of compactness and comfort, allowing you to move freely in any environment.

Comfortable Seating: Take a moment to relax on our comfortable seats, providing you with a convenient spot to rest during your daily activities.

Generous Weight Capacity: Each rollator in our Cirrus Collection is designed to support users weighing up to 136kg, ensuring strength and stability for users of all sizes.

Removable Bag: Fashionable shopping bag with strap.

Stick Holder: Handy holder for a walking stick.



AVAILABLE IN 3 COLOURS

90716 Burnt orange

90717 Dark grey

90718 Purple



CIRRUS+

LIGHTWEIGHT
4 WHEEL ROLLATOR

All the excellent qualities of the Cirrus Rollator 'plus' upgraded features.

Meet the Cirrus Plus Lightweight 4-Wheel Rollator weighing only 6kg. Engineered for effortless mobility, comfort, and confidence. Made from high-quality aluminium, it offers outstanding strength and durability while remaining lightweight and easy to manoeuvre for everyday use.

Enjoy a smooth, stable ride on any surface thanks to the large treaded wheels, featuring upgraded 9" flat-free PU front tyres and durable PP+PVC wheels. Designed for control and stability, the Cirrus Plus helps you move safely and comfortably indoors or outdoors.

Comfort comes standard with a padded seat, height-adjustable padded backrest, and simple one-touch handle height adjustment, allowing you to rest and customise your fit with ease. With a sleek metallic finish, it delivers both modern style and practical performance.

Thoughtfully designed for daily convenience, it includes a large removable zip-closure shopping bag, a walking stick holder, and supports users up to 136kg - making the Cirrus Plus a reliable, stylish mobility solution you can trust.



90743 Dark Grey

Only
6.9kg

Max
Weight
136kg



STRATUS

BARIATRIC ROLLATOR WITH SUSPENSION

Empowering Comfort and Stability for Every Journey!

Especially designed to support users with a max weight of 225kg while providing unparalleled comfort and stability.

Experience the ultimate blend of comfort, stability, and convenience. Don't let mobility challenges hold you back – step into a world of freedom and independence today.

AVAILABLE IN



Pink



Blue



NIMBUS

CARBON FIBRE ROLLATOR

One of the lightest carbon fibre rollators on the market

We are thrilled to introduce the ultimate lightweight Nimbus Carbon Fibre Rollator weighing only 5.2kg.

Crafted with a carbon fibre frame, this rollator combines unparalleled lightness with exceptional robustness, ensuring durability without compromise.

Don't let mobility limitations hold you back – embrace freedom and independence with every step.



90695

GEMINI Fold Down Rollator

The new Gemini Fold Down Rollator – your ultimate companion for mobility and convenience! Say goodbye to bulky mobility aids. Our rollator features a lightweight aluminium frame for easy manoeuvrability without sacrificing durability. Our innovative dual fold mechanism allows the unit to be folded flat effortlessly, making transport and storage a breeze.



Max Weight
136kg

Only
8kg

- Lightweight Aluminium frame
- Innovative dual fold mechanism allows unit to be folded flat for easy transport and storage
- Suitable for indoor and outdoor use
- Simple, one-touch, easily adjustable height adjustable push handles
- Push to lock brakes
- Chunky large treaded castors with Flat Free Tyres
- Removable padded backrest
- Comfortable breathable seat
- Easily removable shopping bag with carry strap
- Gloss black paint finish
- Push handles: 840mm - 980 mm
- Seat height: 620mm. Seat width: 430mm
- L: 720mm W: 610mm
- Folded: 235 x 960 x 400mm
- Max user weight: 136kg

90703

ORION Aluminium Transit Wheelchair with Flip up Arms

Suitable for both indoor and outdoor use, this attendant-propelled wheelchair is an excellent choice for the occasional user.

This chair has been crash tested and can be used for secured transport within a motor vehicle.



- Seat belt included
- Half folding backrest
- Flat Free Tyres
- Compact when folded (750 x 800 x 260mm)
- Front wheels: 200mm, Rear wheels: 300mm
- Seat W: 450mm, Seat D: 400mm, Seat H: 500mm
- Width between arms: 460mm
- Min footrest length: 420mm
Max footrest length: 490mm
- Length inc footrests: 1020 - 1030mm
- Width when open: 560mm
- Width when folded: 260mm
- Height: 1000mm
- Weight: 15kg
- Max user weight: 120kg
- BS EN ISO 10993-10

90681

APOLLO

LIGHTWEIGHT WHEELCHAIR

Available in attendant propelled (AP) or self-propelled (SP) models

- Lightweight Aluminium frame weighing 9.2kg
- Half folding backrest and simple frame folding mechanism. Quick to fold and store in the back of most cars
- Swing away or easily removable leg rests with height adjustable footplates
- The double cross brace gives superb durability and strength
- Comfort is ensured with 18" seat width, PU padded armrests and soft breathable padded seat and backrest upholstery
- Push to lock attendant brakes give the carer full control, easily braking the chair without having to bend down
- Flat free tyres and front castors remove the worry of punctures
- Lap belt and adjustable anti tip bars.

- Width between arms: 45cm
- Seat W:45cm, D:42.5cm, H:45cm
- Armrest height from seat: 21.5cm
- Folded width: 28cm
- Backrest height: 460mm
- Total weight: 9.2kg
- Weight without footrests: 8kg



Attendant Propelled



Self Propelled

mm L: 910
W: 640
H: 910

Only 9.2kg

Max Weight 136kg

| AP | SP |
|-------|-------|
| 90737 | 90747 |

YO-YO

FULL SUSPENSION WHEELCHAIR

Available in attendant propelled (AP) or self-propelled (SP) models

Lightweight and easy to manoeuvre, this folding wheelchair is designed to deliver exceptional comfort, reliability, and convenience for everyday use.


Featuring a full 4-wheel suspension system, it provides a smooth, stable ride across a variety of surfaces, reducing vibration and improving overall user comfort. The generously padded upholstery offers excellent support for longer periods of sitting, while the push-to-lock attendant brakes, adjustable armrests and swing-away leg rests allow for effortless transfers and personalised positioning.

Flat-free tyres ensure peace of mind with zero maintenance and dependable performance wherever you go.

Compact when folded for transport or storage, this wheelchair is ideal for both home and travel.

Available in both attendant and self-propelled models, it offers flexibility to suit individual mobility needs while combining practicality with modern design.

- ✓ 4 Wheel Suspension
- ✓ Lightweight
- ✓ Compact Folding
- ✓ Flat Free Tyres



| AP | SP |
|-------|-------|
| 90741 | 90746 |

AVAILABLE
IN 3 COLOURS

ALTO

LIGHTWEIGHT 3 WHEEL ROLLATOR



- 90750 Purple
- 90751 Orange
- 90752 Grey

Lightweight and easy to manoeuvre, this modern 3-wheel rollator offers confidence and freedom both indoors and outdoors. Its slim design moves smoothly through tight spaces, while the removable storage bag keeps everyday essentials close at hand.

Designed with comfort and convenience in mind, the ergonomic height-adjustable handles ensure a personalised fit, and the self-locking easy-fold mechanism allows for quick, secure folding when transporting or storing the rollator.

Available in three vibrant metallic colours, this rollator combines dependable support with a sleek, contemporary look—perfect for those who want practical mobility assistance without sacrificing style.



3 Wheel Aluminium Rollator

Designed for effortless mobility and everyday confidence. Built with a lightweight yet durable grey powder-coated aluminium frame, this rollator is easy to manoeuvre and simple to fold for transport or storage. Height-adjustable push handles ensure a comfortable, customised fit, while smooth-rolling 200mm PU castor wheels provide stability and control.

Safety and convenience come as standard, with easy push-to-lock brakes and a practical storage bag. Compact yet sturdy, it offers a maximum user weight of 125kg with an adjustable handle height of 800–980mm, this rollator is the perfect blend of comfort, strength, and modern practicality.

- Height: 800 - 980mm
- Length: 670mm
- Total width at rear: 685mm
- Max weight: 125kg
- BS EN ISO 10993-10

90386

LUNA

LIGHTWEIGHT 4 WHEEL ROLLATOR

Lightweight 4-wheel rollator designed for ease of mobility, featuring a three-position adjustable seat height, a comfortably padded seat and supportive backrest, ergonomically adjustable handles for personalized comfort, durable flat-free wheels for smooth and worry-free movement, and a removable shopping bag to carry essentials conveniently.



90758



4 Wheel Aluminium Rollator

This rollator provides stability, comfort, and confidence with every step. Its lightweight aluminium frame is easy to manoeuvre, while the foldable design makes transport and storage simple and hassle-free.

Built for everyday comfort, it features a padded seat and supportive backrest, allowing users to rest whenever needed. The 200mm rubber castors provide smooth rolling and grip across a variety of indoor and outdoor surfaces.

Practical and adjustable, the rollator offers an adjustable handle height from 800–980mm to suit different users. With a strong 145kg weight capacity and a compact footprint, it's a dependable mobility aid for daily use.

- Height: 800 - 980mm
- Length: 670mm
- Width: 630mm
- Width between arms: 445mm
- Seat height: 550mm
- Max weight: 145kg
- BS EN ISO 10993-10



90387

Moving & Handling

Discover our 'Moving and Handling' section, a dedicated collection of products designed to assist occupational therapists and caregivers in providing safe and efficient patient care.

Within this section, you'll find a comprehensive range of equipment and aids specially crafted to facilitate the seamless movement of patients from one place to another. Our selection includes items that prioritize patient comfort, safety, and ease of use, enabling professionals to deliver top-quality care.

With these valuable tools at your disposal, you can confidently and effectively manage patient transfers and handling, ensuring their well-being and comfort at every step.

RS1 Patient Turner

- Enables sit-to stand for semi-ambulant users
- Height adjustable handle
- Height adjustable knee pads
- Slip resistant turning disc
- Easy to use and manoeuvre
- Length: 472mm
- Width: 475mm
- Height: 914mm
- Castor Size: 75mm
- Max weight: 150kg
- BS EN ISO 10993-10



90330



RS2 Patient Transfer Assist

- Low step on height base for ease with clients who have minimal leg movement
- Anti slip base ensures the client is safe and secure in transfer
- Low friction castor for ease of manoeuvrability
- One foot touch brake keeps the base stable when embarking and disembarking from the raiser
- Versatile handle design ensures carers of all shapes and sizes can easily manoeuvre it
- Adjustable knee pad gives the client support and comfort when in use
- Length: 635mm, Width: 607mm, Height: 1186mm
- Max weight: 160kg
- BS EN ISO 10993-10

90329



Ladder Belt (Sliding Support Strap)

- For use on RS2 Patient Transfer Assist
- Length: 2000mm
- Width: 100mm
- Height: 20mm

90359



Patient Handling Belt

- For use on RS2 Patient Transfer
- 3 sizes available: Small, Medium and Large
- Length (mm): S:1050 / M:1300 / L:1500
- Width (mm): S:975 / M:1150 / L:1350
- Height: 100mm

| | |
|--------|--------|
| Small | 90333S |
| Medium | 90333M |
| Large | 90333L |



RS4 Patient Transfer Assist

- Light Weight, Strong and Robust frame
- Easy flip down support pads
- Easy to open leg mechanism. Works in conjunction with most types of seating and chairs
- Enables users to participate in transfers, maintaining patient independence
- Ideal for those users who can raise from a seated position to a standing position
- Adjustable kneepad
- Compact narrow frame. Manoeuvre in small spaces
- Length: 820mm, Width: 530-730mm, Height: 1060mm
- Seat height: 620-670mm
- Max weight: 185kg
- BS EN ISO 10993-10



| | |
|------------------|-------|
| 3" Subassembly | 90580 |
| 3" Full Assembly | 90720 |
| 4" Subassembly | 90724 |
| 4" Full Assembly | 90733 |



AP Aluminium Shower Commode Chair

- Removable PU Back and PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 460mm
- Seat height: 510-535mm
- Height: 985-1010mm
- Length: 860mm
- Width: 569mm
- BS EN ISO 10993-10



90209



SP Aluminium Shower Commode Chair

- Removable PU Back and Horseshoe Seat with Aperture Plug
- 2 x Braked Nylon Castors and 2 x 24" MAG Wheels with PU Tyres and Q.R Pins
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 460mm
- Seat height: 510-535mm
- Height: 985-1010mm
- Length: 860mm
- Width: 569mm
- BS EN ISO 10993-10



90349



AP HD Aluminium Shower Commode Chair

- Removable PU Back
PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 560mm
- Seat height: 510-535mm
- Height: 985-1010mm
- Length: 860mm
- Width: 669mm
- BS EN ISO 10993-10



90350

AQUA LIBRA

STAINLESS STEEL SHOWER CHAIR RANGE



CONTACT OUR SALES TEAM ON
01249 760581
sales@urathon.com

Aqua Libra AP

ATTENDANT PROPELLED

Aqua Libra SP

SELF PROPELLED



3 HANDLE OPTIONS AVAILABLE



| | Straight Handle | | Push Handle | | Pram Handle | |
|------------|-----------------|-------------|-------------|-------------|-------------|-------------|
| | AP | SP | AP | SP | AP | SP |
| 17" | 90646-50001 | 90649-50001 | 90646-50010 | 90649-50010 | 90646-50023 | 90649-50023 |
| 19" | 90647-50002 | 90650-50002 | 90647-50011 | 90650-50011 | 90647-50024 | 90650-50024 |
| 21" | 90648-50003 | 90651-50003 | 90648-50012 | 90651-50012 | 90648-50025 | 90651-50025 |



FEATURES AND BENEFITS

- Fully modular construction – ability to adapt the chair to suit any user
- Stainless steel, white powder coated frame – cannot rust and lends itself well to constant recycling if required
- Height adjustment frame as standard, designed to be able to fit over most standard toilets
- Large commode pan capacity
- 17", 19" and 21" frame width options
- Max user weight: 160Kg
- Available in both attendant and self propelled variants
- Wide range of seat variances – ideal for all toileting requirements
- Adjustable / easily removable arms can be supplied with various width options
- Easily removable backrest with 3 handle options available
- Easy to clean seat and backrest materials
- Fully removable swing away footrests with simple length adjustment and nylon footplates.



AP XXL Aluminium Shower Commode Chair


- Removable PU Back, PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 600mm
- Seat height: 510-535mm
- Length: 860mm
- Width: 715mm
- Height: 985-1010mm
- Max weight: 300kg
- BS EN ISO 10993-10



90513

Self Propel Wheel Kit available

90351



AP Aluminium Shower Commode Chair, Tilt in Space

- Lightweight aluminium frame - Easy to manoeuvre and offers long-term rust-free durability.
- Tilt and space chair – ideal for showering people comfortably
- 4 x Large, braked, robust castors for durability and manoeuvrability.
- Polyurethane seat, seat-infill, back and armrests are water-resistant, easy to clean and comfortable.
- Full height adjustable to suit all heights of the users
- Simple swing away and height adjustable footrests
- Easy to clean seat and removable back
- Flip up arm rests – For easy side transfer
- Large capacity commode pan and lid
- BS EN ISO 10993-10



90581



Chrome Plated Glideabout Commode

- 4 x Castors
- Removable, Grey Aperture Seat with Padded Vanity Cover
- Removable Arm Rests, Backrest, Leg Rests
- Potty
- Grey Vinyl Upholstery
- Length: 840mm
- Width: 525mm
- Height: 985mm
- Seat height: 525mm
- Max weight: 165kg

90099



Height Adjustable Chrome Plated Glideabout Commode

- Height Adjustable
- Height: 985 - 1085mm
- Seat height: 525 - 625mm
- Max weight: 165kg

90294



Mild Steel Glideabout Commode

- 4 x Castors
- Removable, Grey Aperture Seat with Padded Vanity Cover
- Removable Arm Rests, Backrest, Leg Rests
- Potty
- Grey Vinyl Upholstery
- Length: 922mm
- Width: 530mm
- Height: 962mm
- Seat height: 528mm
- Max weight: 165kg
- BS EN ISO 10993-10

90517



Height Adjustable Mild Steel Glideabout Commode

- Height Adjustable
- Height: 985-1085mm
- Max weight: 165kg
- BS EN ISO 10993-10

90675



KD Glideabout Commode

- Easy to Dissassemble
- Removable Aperture
- Padded Vinyl Seat
- Removable Back, Arms, Swing Angle Leg Rests
- Length: 902mm
- Width: 520mm
- Height: 908mm
- Seat height: 530mm
- Max weight: 190kg
- BS EN ISO 10993-10

90519

Hospital Portering Chair

- Full PU Seat and back – Durable, robust, easy to clean and reduces any cross-infection risk. Ideal in heavy usage situations
- Strong balanced footplate – Durable with the ability to take the full user weight without tipping
- The footplate is made from HDPE plastic – Robust and ideal for infection control. Competitors tend to use wood which is a cross infection risk
- Stainless steel frame – Ensure product longevity. Allows the product to be used outside without fear of rust. Easy to wipe surface reduces risk of cross infection. Competitors are made from mild steel that can rust when exposed to the inclement weather
- 300kg SWL – Most competitors are only 190 KGs, so this can be used with the heaviest of patients
- Moulded drop-down arms – Allows for quick and safe side transfers
- Accessory mounts – Allows IV poles etc to be used with the chair
- Xray/notes pocket as standard – Enables the users records to be transported safely with the user from department to department.
- BS EN ISO 10993-10

- Size: 600mm x 865mm x 1040mm
- Seat Width: 540mm
- Seat Depth: 410mm
- Seat Height: 540mm
- Backrest Height: 480mm



90337



Aquarius KD Shower Commode Chair

- Can be easily used as a; shower chair, freestanding toilet raiser and commode
- Removable padded backrest
- Removable arms
- Easily removable height adjustable legs
- Non slip rubber feet
- Removable potty
- Aperture seat with padded vanity cover
- Length: 395mm
- Width: 505mm
- Weight: 750 – 877mm
- Seat Height: 420 – 545mm
- Seat width: 460mm
- Weight: 5.1kg
- Max user weight: 130kg

90667



Greenwall KD Ali Shower Chair

- Removable, padded armrests for comfort and added security
- Blow-molded seat and back provide comfort and strength and is warm to the touch
- Drainage holes in bench reduce slipping
- Aluminium frame
- Leg height adjusts in 25mm increments
- Product Specifications
- Product Weight: 4.4kg
- Max user weight: 136kg
- Seat width: 400mm
- Seat width : 330mm
- Seat to Floor Height: 400mm – 530mm
- Width Between Arms: 490mm

90666



Shower Stool

- Height Adjustable
- Powder Coated Steel Frame
- Easy to remove PP Seat
- Aluminium option
- Depth: 370-395mm
- Width: 495-540mm
- Seat width: 380mm
- Seat depth: 280mm
- Seat height: 480-630mm
- Width between arms: 485mm
- Max weight: 160kg
- BS EN ISO 10993-10

90115



Drop Down Rail with Support Leg

- Length: 692mm
- Width: 100mm
- Height: 750mm
- Max weight: 190kg
- BS EN ISO 10993-10

90146



HD Suction Grab Rail 44/56cm

- Bar adjustable over 12cm
- Strong rubber suction cups
- Secure suction cups by lowering the lever
- Max weight: 40 kg

90631



HD Suction Grab Rail 55/68cm

- Bar adjustable over 12cm
- Strong rubber suction cups
- Secure suction cups by lowering the lever
- Max weight: 40 kg

90632



HD Suction Grab Rail 65/78cm

- Bar adjustable over 12cm
- Strong rubber suction cups
- Secure suction cups by lowering the lever
- Max weight: 40 kg

90633



2"

Raised Toilet Seat, No Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90635



4"

Raised Toilet Seat, No Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90636



6"

Raised Toilet Seat, No Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90637



2"

Raised Toilet Seat, With Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90660



4"

Raised Toilet Seat, With Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90661



6"

Raised Toilet Seat, With Lid

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90662



Toilet Surround Height and Width Adjustable

- Height & Width Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 584 - 734mm
- Height: 645 - 795mm
- Max weight: 190kg
- BS EN ISO 10993-10

90438



Toilet Surround

- Height Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 594mm
- Height: 645 - 795mm
- Max weight: 180kg
- BS EN ISO 10993-10

90702



Static Shower Chair

- Height Adjustable
- Length: 480-505mm
- Width: 620 - 645mm
- Height: 775-925mm
- Seat height: 440-590mm
- Max weight: 160kg
- BS EN ISO 10993-10

90300



4 Part Bath Step

- Length: 375mm
- Width: 480mm
- Height: 110mm
- Max weight: 250kg

90116



Male Urinal Bottle with lid

- Length: 126mm
- Width: 98mm
- Height: 225mm

90597

90040
6"

90041
8"

90042
12"



Bath Seat

- Strong suction feet
- Available in 150mm, 250mm and 350mm height options
- Length: 290mm
- Width: 450mm
- Max weight: 190kg

90398



Swivel Bath Stool

- PU Seat with Lockable Swivel Function
- Storage Tray
- Height Adjustable
- Length: 380mm
- Width: 380mm
- Height: 410 - 580mm
- Max weight: 135kg



Toilet Frame, Height Adjustable

- Injection Moulded Seat
- L: 443-471mm
- W: 607-612mm
- H: 622-771mm
- Seat H: 428-578mm
- Max weight: 190kg
- BS EN ISO 10993-10

90341



Toilet Frame, Width & Height Adjustable

- Injection Moulded Seat
- L: 443-471mm
- W: 570-720mm
- H: 622-771mm
- Seat H: 428-578mm
- Max weight: 190kg
- BS EN ISO 10993-10

90532



Swivelling Bath Chair

- Adjustable Width
- Width: 770 - 820mm
- Height: 530mm
- Width between arms: 600-700mm
- Max weight: 150kg
- BS EN ISO 10993-10

90071



Toilet Frame, Clip on Seat & Chute

- Height Adjustable
- Length: 465 - 495mm
- Width: 558mm
- Height: 690 - 835mm
- Seat height: 483-628mm
- Max weight: 190kg
- BS EN ISO 10993-10

Version with lid available (90698)

90120



26" Bath Board 4 Slat

- Adjustable securing brackets
- Length: 240mm
- Width: 600mm
- Height: 20mm

90210



27" Bath Board 4 Slat

- Adjustable securing brackets
- Length: 240mm
- Width: 690mm
- Height: 20mm

90211



28" Bath Board 4 Slat

- Adjustable securing brackets
- Length: 240mm
- Width: 710mm
- Height: 20mm

90212



26" Bath Board 6 Slat

- Adjustable securing brackets
- Length: 360mm
- Width: 660mm
- Height: 20mm
- Max weight: 190kg

90254



27" Bath Board 6 Slat

- Adjustable securing brackets
- Length: 360mm
- Width: 690mm
- Height: 20mm
- Max weight: 190kg

90255



28" Bath Board 6 Slat

- Adjustable securing brackets
- Length: 360mm
- Width: 710mm
- Height: 20mm
- Max weight: 190kg

90256

yuwell Anytime™

CONTINUOUS GLUCOSE MONITORING SYSTEM

A world-first rechargeable CGM transmitter

With readings every three minutes and unlimited data storage, this cutting-edge device uses Bluetooth technology to sync seamlessly with your phone or tablet. The Anytime CT-3 CGM transmitter lasts an impressive two years, reducing both costs and environmental waste. Its state-of-the-art ergonomic design ensures user-friendly operation, while superior sterile packaging enhances safety. Best of all, you only need to replace the sensor every two weeks - keeping the transmitter for long-term use!



NO FINGER
PRICKS



NO SCANNING



COMFORTABLE
TO WEAR FOR
14 DAYS



CUSTOMISABLE
ALERTS & ALARMS



RECHARGEABLE
TRANSMITTER
LASTS 2 YEARS

Manage your glucose Anytime, Anywhere!



MARD
6.9%



Order Now
yuwellanytime.co.uk

Bedroom Aids

Welcome to our 'Bedroom Aids' section, where you'll find a curated selection of essential products designed to enhance safety, comfort, and convenience in your bedroom. Explore our range of bed rails, bed supports, and commodes, thoughtfully chosen to provide the support you need to create a secure and comfortable bedroom environment. Whether you're looking for assistance with mobility or seeking solutions to improve your overall quality of life, our bedroom aids are here to help you rest easy.



Bedstick Slatted Bed

- Length: 822mm
- Width: 622mm
- Height: 500mm
- Max weight: 190kg
- BS EN ISO 10993-10

90321



Double Bedstick

- Adjustable Width
- Powder Coated White
- Length: 335mm
- Width: 960 - 1555mm
- Height: 440mm
- Max weight: 140kg
- BS EN ISO 10993-10

90354



Double Bedstick

- Adjustable Width
- Chrome Plated
- Length: 335mm
- Width: 960 - 1555mm
- Height: 440mm
- Max weight: 140kg

90355



Bedstick For Divan Bed

- Length: 805mm
- Width: 622mm
- Height: 646mm
- Max weight: 190kg
- BS EN ISO 10993-10

90322



Straps (PAIR) for UEL90321/90322

- Length: 2500mm
- Width: 25mm
- Max weight: 50kg

90324



Bed Grab Rail

- Length: 653mm
- Width: 435mm
- Height: 510mm
- Max weight: 165kg
- BS EN ISO 10993-10

90521



Bed Lever with Straps

- Length: 610mm
- Width: 362mm
- Height: 561mm
- Max weight: 160kg
- BS EN ISO 10993-10

90325



Stacking Commode Chair

- Height Adjustable
- Fixed Arms
- HF Vinyl Welded Padded Seat & Back
- Brown Powder Coated Steel Frame
- Anti Rattle Ring & DIP Moulded Ferrule
- Commode Pan & Lid
- Length: 525-545mm
- Width: 551mm
- Height: 785-930mm
- Seat Height: 452-602mm
- Max weight: 160kg
- BS EN ISO 10993-10

90096



Stacking Commode Chair with Removable Arms

- Height Adjustable
- Removable Arms
- HF Welded Padded Seat & Back
- Brown Padding and Brown Powder Coated Steel
- Anti Rattle Ring & DIP Moulded Ferrule
- Commode Pan & Lid
- Length: 525-545mm
- Width: 551mm
- Height: 785-930mm
- Seat Height: 452-602mm
- Max weight: 160kg
- BS EN ISO 10993-10

90323



Bed Bumpers with Mesh Divide

- Blue with black mesh
- Length: 2000mm
- Width: 710mm
- Height: 10mm

90522



Bed Bumpers with Mesh Divide

- Cream with white mesh
- Length: 2000mm
- Width: 710mm
- Height: 10mm

90560



Postural Wedge

- Non slip material
- Length: 530mm
- Width: 275mm
- Height: 210mm

90397

Freedom 2 Bed

- 4 section profiling bed
- Electrically operated backrest, height and knee break angle
- See Freedom 2 manual for details
- Available with and without side rails
- Length: 2230mm
- Width: 980mm
- Height: 335 - 735mm
- Max weight: 178kg
- BS EN ISO 10993-10



90533



Freedom 2 Bed Lever

- Compatible with the Freedom 2 Bed
- Length: 126mm
- Width: 286mm
- Height: 713mm
- Max weight: 160kg

90558

Around The Home

In our 'Around the Home' range you'll find a thoughtfully curated collection of products designed to improve the functionality and comfort of your living space. Explore our range of footstools, overbed tables, and perching stools, each carefully selected to enhance your daily life. Whether you're seeking extra support while sitting or assistance with daily tasks, our home-focused aids are here to make your surroundings more accommodating and user-friendly. Elevate your home environment and enjoy increased convenience and ease with our selection of purposeful products.



- Height Adjustable
- Width: 1010mm
- Depth: 480mm
- Height: 710-910mm
- Max Weight: 5kg

Overchair Table

This curved overchair/bed table is both width and height adjustable allowing the table to be set up perfectly to the requirements of the user. This table is great for use with chairs, single beds and wheelchairs.

The kidney shaped top is designed to curve gently around the user reducing the need for them to reach forward, it also allows them to rest their arms comfortably on the table top whilst still having a large usable work area.

90615



Overbed Table

- Height Adjustable
- NO Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Top
- Length: 697mm
- Width: 413mm
- Height: 630 - 930mm
- BS EN ISO 10993-10

90290



Overbed Table with Castors

- Height Adjustable
- 4 x 2" Braked Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Table Top
- L: 697mm. W: 413mm
- Height: 700 - 1000mm
- BS EN ISO 10993-10

90291



Vinyl Footstool

- Height Adjustable
- 4 x Ferrule
- Length: 500mm
- Width: 315mm
- Height: 380-530mm
- Seat height: 380-535mm
- Max weight: 100kg
- BS EN ISO 10993-10

90454



Vinyl Footstool on Castors

- Height Adjustable
- 4 x 40mm Castors
- Length: 500mm
- Width: 345mm
- Height: 375-525mm
- Seat height: 375-525mm
- Max weight: 100kg
- BS EN ISO 10993-10

90455

Crutch Caddy

GIVING INDEPENDENCE

2
FREE
RUBBER
CRUTCH
TIPS

An innovative mobility aid to safely store a pair of crutches in an upright position ready to use.



Easy to use
lock and release



Suits all
standard crutches*



Glow in the
dark

Home

Workplace

Hospitals

Waiting Rooms

Functions

Keeps crutches
tidy and accessible



**2
FREE
RUBBER
CANE
TIPS**

Cane Caddy

GIVING INDEPENDENCE



An innovative mobility aid to safely store a cane in an upright position ready to use.



Available in two hole sizes*



Lightweight and portable

The unique Crutch and Cane Caddy may be positioned wherever the safe storage of crutches or canes are needed. By your bedside, in the bathroom or next to a comfortable chair. Lightweight and portable, they provide easy access to crutches or canes, giving independence and freeing up carer's time.

Keeps canes tidy and accessible



Perching Stool

- HF Vinyl Padded Seat
- Adjustable Height
- Length: 415- 460mm
- Width: 561 - 620mm
- Height: 561 - 705mm
- Seat height: 540-686mm
- Max weight: 160kg
- BS EN ISO 10993-10

90167



Perching Stool with Arms

- HF Padded Seat
- Steel Loop Arms, No Back
- Adjustable Height
- Length: 415- 460mm
- Width: 561 - 620mm
- Height: 561 - 705mm
- Seat height: 540-686mm
- Max weight: 160kg
- BS EN ISO 10993-10

90168



Perching Stool with Arms and Back

- HF Vinyl Seat
- Steel Loop Arms & Steel Back
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 815 - 960mm
- Seat height: 540 - 685mm
- Max weight: 160kg
- BS EN ISO 10993-10

90169



Perching Stool with Arms and HF Back

- HF Vinyl Seat
- Steel Loop Arms & HF Back
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 830 - 975mm
- Seat height: 540 - 685mm
- Max weight: 160kg
- BS EN ISO 10993-10

90245



Bariatric Perching Stool with Arms & HF Back

- HF Seat and HF Back
- Adjustable Height
- Length: 500 - 550mm
- Width: 690 - 750mm
- Height: 950 - 1095mm
- Seat height: 570-715mm
- Max weight: 290kg
- BS EN ISO 10993-10

90327



Perching Stool with PU Seat

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 545 - 685mm
- Seat height: 535 - 680mm
- Max weight: 160kg
- BS EN ISO 10993-10

90270



Perching Stool with PU Seat and Steel Arms

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 723-868mm
- Seat height: 535 - 680mm
- Max weight: 160kg
- BS EN ISO 10993-10

90271



Perching Stool with PU Seat, Steel Arms and Back

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 815 - 960mm
- Seat height: 535 - 680mm
- Seat width: 325mm
- Seat depth: 265mm
- Max weight: 160kg
- BS EN ISO 10993-10

90272



Perching Stool with PU Seat, Steel Arms and Padded Back

- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 855 - 1000mm
- Seat height: 535 - 680mm
- Seat width: 325mm
- Seat depth: 265mm
- Max weight: 160kg
- BS EN ISO 10993-10

90273



Perching Stool with PU Seat, Steel Arms and Padded Back

- 35mm Extended Seat
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 855 - 1000mm
- Seat height: 535 - 680mm
- Seat width: 368mm
- Seat depth: 275mm
- Max weight: 160kg
- BS EN ISO 10993-10

90674



Bariatric Perching Stool with PU Seat, PU Back & Arms

- Adjustable Height
- Length: 465-515mm
- Width: 650-710mm
- Height: 790 - 935mm
- Seat height: 485-630mm
- Max weight: 290kg
- BS EN ISO 10993-10

90346



KD Perching Stool

- Removable Arms, Legs and Back
- Steel Arms
- Height Adjustable
- Max weight: 160kg
- BS EN ISO 10993-10

90626



Two Tier Trolley with Brakes

- Push to lock brakes
- Length: 530mm
- Width: 595mm
- Height: 770 - 910mm
- Castor size: 150mm
- Max weight: 125kg
- BS EN ISO 10993-10

90515



2 Tier KD Dinner Trolley

- Height Adjustable
- 2 x Removable Shelves
- Length: 501mm
- Width: 445mm
- Height: 845 - 945mm
- BS EN ISO 10993-10

90072



Half Step

- Adjustable Feet
- Removable slip resistant mat
- Plastic - will not rust
- Length: 430mm, Width: 735mm, Height: 112mm
- Max weight: 190kg

90141



Back Rest

- Folds flat for easy storage
- Sturdy frame with breathable mesh fabric
- Adjustable soft pillow
- 5 levels of angle adjustment
- Width: 580mm, Depth: 515mm, Height: 560mm
- Max weight: 127kg

90668



Urathon[®]

UK WAREHOUSING

Our dedicated UK warehouse allows for quick dispatch to you. No overseas shipping delays. Receive your products quickly, on time and when you need them. That's our promise to you.

SUPPLY FROM STOCK

We fulfil product orders quickly from stock within our self managed UK warehouse. This allows us to always provide fast and reliable turnaround times with next day delivery available.

OUTSTANDING CUSTOMER CARE

Customer care is at the heart of our business. Your satisfaction is paramount and we always aim to go above and beyond.

Call **01249 760581**

Email **sales@urathon.com**

Visit **www.urathon.com**

Urathon[®]

**SPECIALIST SUPPLIERS TO THE MOBILITY
& REHABILITATION INDUSTRY**