

Telephone **01249 760581** or email **sales@urathon.com**

2026

# Product Catalogue



# Urathon<sup>®</sup>

SPECIALIST SUPPLIERS TO THE MOBILITY  
& REHABILITATION INDUSTRY

# Urathon®

Urathon have been designing and brand manufacturing for many of the healthcare industries leading suppliers for almost 30 years including many of the UKs largest providers to Nursing homes, Hospitals, NHS community equipment stores and Mobility retailers.

We offer a complete seamless supply chain solution for all your Healthcare and Mobility product requirements either from our large fully stocked warehousing in Wiltshire, UK or by full mixed containers under your own brand directly from our wholly owned manufacturing base in China.



# CONTENTS

# THOR

## ULTRA - LIGHT FOLDING SCOOTER


Discover effortless mobility with this ultra-light folding scooter, combining sleek modern design with a full carbon-fibre easy-fold frame weighing just 13kg. Exceptionally portable and simple to store, it's ideal for everyday use and travel.

Comfort features include padded upholstery and fold-up, removable armrests, while clear, easy-to-use controls offer horn, speed adjustment, and indicators for confident riding. Quiet 180W motors provide smooth performance and a range of up to 15.5 miles.

Two lightweight, removable aircraft-safe lithium batteries, USB-C charging, and an optional magnetic phone holder complete this smart, convenient mobility solution designed to keep you moving with ease.



**5 YEARS WARRANTY**  
Terms and Conditions apply



Lightweight Carbon Fibre Finish

Only **13kg**

Max Weight **150kg**

+ Range **15.5 miles**

90760

# HALO

## FLAT-FOLDING SCOOTER



Terms and Conditions apply

### Freedom, Comfort, and Convenience

Experience a new level of independence with this stylish and practical folding mobility scooter, designed to make everyday travel easier and more enjoyable. With its modern, sleek design and lightweight construction, it combines portability, comfort, and smart technology in one compact package.

### Designed for Easy Travel

The scooter features a clever fold-flat design, allowing it to be stored quickly and conveniently. Weighing just 16.5 kg including the battery, it can be easily lifted and loaded into a car, making it perfect for days out, shopping trips, or holidays.

### Reliable Performance


Powered by a high-quality aircraft-safe battery, the scooter offers a range of up to 11 miles on a single charge, providing plenty of distance for everyday journeys. For even greater flexibility, the scooter has the option to use two batteries, allowing you to extend your travel range when needed.

### Smart and Simple to Use

Equipped with NFC technology, the scooter offers modern convenience and added security. The simple, easy-to-use controls ensure smooth operation, making it suitable for both new and experienced users.

### Comfort Built In

Featuring adjustable handlebars to suit your preferred driving position, along with a soft upholstered seat and backrest that provide excellent support during longer journeys.



Stylish Satin  
Carbon Fibre Finish

90762

# HERO

## CARBON FIBRE POWERCHAIR

Revolutionary frame design manufactured from aviation grade lightweight carbon fibre, available in stunning ruby red and matt onyx black.

Easy to fold and store, luxury padded backrest and memory foam cushion.

Advanced, intuitive easy to use controls, fold up armrests for easy access.

2x 250w powerful brushless motors powered by an easy removable aircraft-safe lithium battery gives an impressive 10 mile range, capable of climbing 12 degree slopes with 7cm ground clearance.



Terms and Conditions apply



Matt Carbon Fibre Finish



# ZEUS


## AUTO-FOLDING POWERCHAIR



The Zeus Auto Folding Powerchair is a stylish, travel-ready mobility solution built for ultimate convenience. Its standout auto-folding mechanism allows the chair to fold in seconds at the touch of a button, making storage and transport effortless. Once folded, the retractable carry handle ensures easy lifting and smooth handling, whether placing it in a car boot or storing it neatly at home.

Crafted with a lightweight aluminium frame and finished in a sleek Satin Carbon Fibre finish, The Zeus combines durability with modern design. A removable, aircraft-safe lithium battery delivers up to a 6 mile range for confident everyday travel. Front suspension provides a smoother ride, while padded upholstery enhances comfort and support.

With a maximum user weight of 120kg,  
The Zeus offers dependable performance  
in a compact, easy-to-transport package.



Stylish Satin  
Carbon Fibre Finish



90740

# NIRA


## LIGHTWEIGHT POWERCHAIR

This lightweight powerchair delivers excellent comfort, stability, and everyday practicality. Featuring a quality upholstered padded seat and backrest, front suspension, and ruggedised tyres, it's designed to perform smoothly both indoors and outdoors.

With dual removable batteries providing a 12–13 mile range, this powerchair folds easily for transport and storage and can be quickly lifted into most vehicles. Flip-up armrests, swing-away leg supports, anti-tip bars, and airline-friendly batteries make it a reliable and versatile all-round mobility solution.



Terms and Conditions apply



- ✔ Front Suspension
- ✔ Anti-tip and Lap Belt
- ✔ Compact Folding
- ✔ Ruggedised Tyres



90759

# AURA


## ULTRA LIGHTWEIGHT POWERCHAIR



Terms and Conditions apply



- ✓ Ultra Lightweight
- ✓ Anti-tip and Lap Belt
- ✓ Compact Folding
- ✓ Bluetooth Remote
- ✓ Aircraft-safe Battery



Stylish Satin  
Carbon Fibre Finish



Experience effortless mobility with this ultra-lightweight powerchair, weighing just 12kg. Combining stability, comfort, and contemporary design, it's perfectly suited for everyday use. The chair folds quickly and easily, making it ideal for lifting into most vehicles and convenient for compact storage at home or on the go.

Featuring advanced Bluetooth remote control, it offers seamless operation without the need for wiring or programming. Simply connect and go for smooth, responsive handling.

Designed with comfort in mind, it includes a generously padded seat and backrest, flip-up armrests for easy transfers, and adjustable leg plates for personalised support.

With an aircraft-safe battery, an impressive range of up to 10 miles, and built-in anti-tip bars, this powerchair delivers confidence, safety, and convenience wherever your journey takes you.

90755



### LEO Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Breathable fabric seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large 10" Front and 8" Rear wheels with PU Flat Free Tyres
- High quality easily removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Max Weight: 130kg

Small	90653
Medium	90652



### ARIES Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Uprated soft padded seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large uprated 10" Front and 8" Rear wheels with PU Flat Free Tyres
- High quality easily removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Seat Height: Small 540mm, Medium 620mm
- Max Weight: 130kg

Small	90655
Medium	90654



### TAURUS Rollator

- One touch adjustable height push handles with Ergonomic comfort grips
- Luxury Genuine stitched leather finishing on padded seat
- Height and depth adjustable padded backrest
- Push to lock brakes with safe internal routed brake cables
- Large uprated 10" Front and 8" Rear wheels with chunky treaded PU Flat Free Tyres
- Genuine stitched leather removable zip-shut storage bag
- Easy folding with simple click to lock mechanism
- Simple to attach cane holder
- Weight: Small 7.5kg, Medium 7.8kg
- Length: 720mm
- Width: Folded 280mm, Unfolded 610mm
- Seat Height: Small 540mm, Medium 620mm
- Max Weight: 130kg



Stitched Leather

Small	90657
Medium	90656

# CIRRUS

## LIGHTWEIGHT 4 WHEEL ROLLATOR

Discover the perfect blend of elegance and practicality with our Cirrus Collection of foldable aluminium rollators.

**Lightweight Aluminium Frame:** Engineered from high-quality aluminium, our rollators offer exceptional strength and durability while remaining lightweight for easy manoeuvrability.

**Smooth Rolling Wheels:** Glide effortlessly on our PP+PVC 8" wheels, designed to provide a smooth and stable ride on any surface.

**Sleek Metallic Finish:** Make a statement with our rollators' sleek metallic finish, combining sophistication with practicality.

**Compact Design:** With dimensions tailored for easy navigation, our rollators offer the perfect balance of compactness and comfort, allowing you to move freely in any environment.

**Comfortable Seating:** Take a moment to relax on our comfortable seats, providing you with a convenient spot to rest during your daily activities.

**Generous Weight Capacity:** Each rollator in our Cirrus Collection is designed to support users weighing up to 136kg, ensuring strength and stability for users of all sizes.

**Removable Bag:** Fashionable shopping bag with strap.

**Stick Holder:** Handy holder for a walking stick.



### AVAILABLE IN 3 COLOURS

90716 Burnt orange

90717 Dark grey

90718 Purple



# CIRRUS+

## LIGHTWEIGHT 4 WHEEL ROLLATOR



90743 Dark Grey

Only  
6.9kg

Max  
Weight  
136kg


### All the excellent qualities of the Cirrus Rollator 'plus' upgraded features.

Meet the Cirrus Plus Lightweight 4-Wheel Rollator weighing only 6kg. Engineered for effortless mobility, comfort, and confidence. Made from high-quality aluminium, it offers outstanding strength and durability while remaining lightweight and easy to manoeuvre for everyday use.

Enjoy a smooth, stable ride on any surface thanks to the large treaded wheels, featuring upgraded 9" flat-free PU front tyres and durable PP+PVC wheels. Designed for control and stability, the Cirrus Plus helps you move safely and comfortably indoors or outdoors.

Comfort comes standard with a padded seat, height-adjustable padded backrest, and simple one-touch handle height adjustment, allowing you to rest and customise your fit with ease. With a sleek metallic finish, it delivers both modern style and practical performance.

Thoughtfully designed for daily convenience, it includes a large removable zip-closure shopping bag, a walking stick holder, and supports users up to 136kg - making the Cirrus Plus a reliable, stylish mobility solution you can trust.



# STRATUS

## BARIATRIC ROLLATOR WITH SUSPENSION

**Empowering Comfort and Stability for Every Journey!**

Especially designed to support users with a max weight of 225kg while providing unparalleled comfort and stability.

Experience the ultimate blend of comfort, stability, and convenience. Don't let mobility challenges hold you back – step into a world of freedom and independence today.

**AVAILABLE IN**

90734 Pink

90735 Blue

Max Weight  
225kg

24"  
Seat Width



# NIMBUS

## CARBON FIBRE ROLLATOR

**One of the lightest carbon fibre rollators on the market**

We are thrilled to introduce the ultimate lightweight Nimbus Carbon Fibre Rollator weighing only 5.2kg.

Crafted with a carbon fibre frame, this rollator combines unparalleled lightness with exceptional robustness, ensuring durability without compromise.

Don't let mobility limitations hold you back – embrace freedom and independence with every step.



90695

## GEMINI Fold Down Rollator

The new Gemini Fold Down Rollator – your ultimate companion for mobility and convenience! Say goodbye to bulky mobility aids. Our rollator features a lightweight aluminium frame for easy manoeuvrability without sacrificing durability. Our innovative dual fold mechanism allows the unit to be folded flat effortlessly, making transport and storage a breeze.



Max Weight  
136kg

Only  
8kg

- Lightweight Aluminium frame
- Innovative dual fold mechanism allows unit to be folded flat for easy transport and storage
- Suitable for indoor and outdoor use
- Simple, one-touch, easily adjustable height adjustable push handles
- Push to lock brakes
- Chunky large treaded castors with Flat Free Tyres
- Removable padded backrest
- Comfortable breathable seat
- Easily removable shopping bag with carry strap
- Gloss black paint finish
- Push handles: 840mm - 980 mm
- Seat height: 620mm. Seat width: 430mm
- L: 720mm W: 610mm
- Folded: 235 x 960 x 400mm
- Max user weight: 136kg

90703

## ORION Aluminium Transit Wheelchair with Flip up Arms

Suitable for both indoor and outdoor use, this attendant-propelled wheelchair is an excellent choice for the occasional user.

This chair has been crash tested and can be used for secured transport within a motor vehicle.



- Seat belt included
- Half folding backrest
- Flat Free Tyres
- Compact when folded (750 x 800 x 260mm)
- Front wheels: 200mm, Rear wheels: 300mm
- Seat W: 450mm, Seat D: 400mm, Seat H: 500mm
- Width between arms: 460mm
- Min footrest length: 420mm  
Max footrest length: 490mm
- Length inc footrests: 1020 - 1030mm
- Width when open: 560mm
- Width when folded: 260mm
- Height: 1000mm
- Weight: 15kg
- Max user weight: 120kg
- BS EN ISO 10993-10

90681

# APOLLO

## LIGHTWEIGHT WHEELCHAIR

Available in attendant propelled (AP) or self-propelled (SP) models

- Lightweight Aluminium frame weighing 9.2kg
- Half folding backrest and simple frame folding mechanism. Quick to fold and store in the back of most cars
- Swing away or easily removable leg rests with height adjustable footplates
- The double cross brace gives superb durability and strength
- Comfort is ensured with 18" seat width, PU padded armrests and soft breathable padded seat and backrest upholstery
- Push to lock attendant brakes give the carer full control, easily braking the chair without having to bend down
- Flat free tyres and front castors remove the worry of punctures
- Lap belt and adjustable anti tip bars.

- Width between arms: 45cm
- Seat W:45cm, D:42.5cm, H:45cm
- Armrest height from seat: 21.5cm
- Folded width: 28cm
- Backrest height: 460mm
- Total weight: 9.2kg
- Weight without footrests: 8kg



Attendant Propelled



Self Propelled



**mm** L: 910  
W: 640  
H: 910

Only 9.2kg

Max Weight 136kg

AP	SP
90737	90747

# YO-YO

## FULL SUSPENSION WHEELCHAIR

Available in attendant propelled (AP)  
or self-propelled (SP) models

Lightweight and easy to manoeuvre, this folding wheelchair is designed to deliver exceptional comfort, reliability, and convenience for everyday use.



Featuring a full 4-wheel suspension system, it provides a smooth, stable ride across a variety of surfaces, reducing vibration and improving overall user comfort. The generously padded upholstery offers excellent support for longer periods of sitting, while the push-to-lock attendant brakes, adjustable armrests and swing-away leg rests allow for effortless transfers and personalised positioning.

Flat-free tyres ensure peace of mind with zero maintenance and dependable performance wherever you go.

Compact when folded for transport or storage, this wheelchair is ideal for both home and travel.

Available in both attendant and self-propelled models, it offers flexibility to suit individual mobility needs while combining practicality with modern design.

- ✓ 4 Wheel Suspension
- ✓ Lightweight
- ✓ Compact Folding
- ✓ Flat Free Tyres



AP	SP
90741	90746

# ZEPHYR

## CRASH TESTED WHEELCHAIR



Terms and Conditions apply



**Attendant  
Propelled**

Meet the Zephyr Wheelchair — where safety, comfort, and performance come together in one sleek design. Built from lightweight yet durable aluminium and fully crash-tested for peace of mind, the Zephyr delivers strength without unnecessary bulk. Whether you're navigating daily routines

or traveling with confidence, its flat-free tyres ensure smooth, worry-free mobility every time. The padded upholstery and additional seat cushion provide superior comfort, making extended use feel effortless.



**Self Propelled**

Designed with flexibility in mind, the Zephyr is available in both attendant-propelled and self-propelled models. Adjustable armrests allow for a tailored fit, enhancing comfort and accessibility throughout the day. Every detail has been carefully considered to support independence, ease of use, and long-term reliability.

Backed by a 5-year warranty, the Zephyr Wheelchair is an investment in quality you can trust. Choose a wheelchair that doesn't just meet expectations — it exceeds them. Experience confidence, comfort, and control with Zephyr.

AP	SP
90742	90748



### 3 Wheel Aluminium Rollator

Designed for effortless mobility and everyday confidence. Built with a lightweight yet durable grey powder-coated aluminium frame, this rollator is easy to manoeuvre and simple to fold for transport or storage. Height-adjustable push handles ensure a comfortable, customised fit, while smooth-rolling 200mm PU castor wheels provide stability and control.

Safety and convenience come as standard, with easy push-to-lock brakes and a practical storage bag. Compact yet sturdy, it offers a maximum user weight of 125kg with an adjustable handle height of 800–980mm, this rollator is the perfect blend of comfort, strength, and modern practicality.

- Height: 800 - 980mm
- Length: 670mm
- Total width at rear: 685mm
- Max weight: 125kg
- BS EN ISO 10993-10

90386

### 4 Wheel Aluminium Rollator

This rollator provides stability, comfort, and confidence with every step. Its lightweight aluminium frame is easy to manoeuvre, while the foldable design makes transport and storage simple and hassle-free.

Built for everyday comfort, it features a padded seat and supportive backrest, allowing users to rest whenever needed. The 200mm rubber castors provide smooth rolling and grip across a variety of indoor and outdoor surfaces.

Practical and adjustable, the rollator offers an adjustable handle height from 800–980mm to suit different users. With a strong 145kg weight capacity and a compact footprint, it's a dependable mobility aid for daily use.

- Height: 800 - 980mm
- Length: 670mm
- Width: 630mm
- Width between arms: 445mm
- Seat height: 550mm
- Max weight: 145kg
- BS EN ISO 10993-10



90387

Lightweight and easy to manoeuvre, this modern 3-wheel rollator offers confidence and freedom both indoors and outdoors. Its slim design moves smoothly through tight spaces, while the removable storage bag keeps everyday essentials close at hand.

Designed with comfort and convenience in mind, the ergonomic height-adjustable handles ensure a personalised fit, and the self-locking easy-fold mechanism allows for quick, secure folding when transporting or storing the rollator.

Available in three vibrant metallic colours, this rollator combines dependable support with a sleek, contemporary look—perfect for those who want practical mobility assistance without sacrificing style.

# ALTO

## LIGHTWEIGHT 3 WHEEL ROLLATOR

AVAILABLE  
IN 3 COLOURS



- 90750 Purple
- 90751 Orange
- 90752 Grey



# LUNA

## LIGHTWEIGHT 4 WHEEL ROLLATOR

Lightweight 4-wheel rollator designed for ease of mobility, featuring a three-position adjustable seat height, a comfortably padded seat and supportive backrest, ergonomically adjustable handles for personalized comfort, durable flat-free wheels for smooth and worry-free movement, and a removable shopping bag to carry essentials conveniently.



Max  
Weight  
135kg

90758

# STAIRLIFT IN A BOX

Introducing our powered footplate and  
swivel seat stairlift for the retail market

Compact • Comfortable • Dependable

With its sleek, matt silver, anodised-rail and comfortable ergonomic seat, a Urathon stairlift looks as good as it performs. Discreet, compact and contemporary, it blends effortlessly into almost any interior without dominating your hallway or staircase.

An LCD display provides clear, simple operating instructions whilst the auto-recharging, battery-powered motor delivers smooth and stable movement with a reassuringly quiet operation, so it won't disturb your household - even during the night.

Suitable for most standard straight staircases (25° to 55° rail angles) and supporting users up to 120kg, this stairlift offers a dependable, elegant and practical way to stay independent in the home you love - all at the touch of a button.



Max  
Weight  
120kg

## Specification

- Comes ready to install and use
- Certified to EN 81-40 standards
- Supplied in one box (including 2.5m dual rails)
- Installation & User manuals
- Supports persons up to 120kg
- Fits rail angles 25° to 55°
- Compact, fold-away design
- Quiet motor, for use day & night
- Auto-recharging dual 12v batteries (works even during power cuts)
- Retractable self-locking seatbelt
- Powered swivel seat & footplate
- Left- & right-hand configuration

For more information, call us today.

**+44 (0)1249 760 581**

[sales@urathon.com](mailto:sales@urathon.com)

**Urathon®**



## RS1 Patient Turner

The RS1 Patient Turner is expertly designed to make transfers safer, smoother, and more comfortable for individuals who are weight-bearing and have good upper body strength. By allowing users to stand and pivot between seated surfaces, the RS1 promotes greater independence while reducing strain for both the user and the caregiver. Simply pulling on the sturdy handlebar helps the user rise into a standing position, enabling a controlled and secure turning movement between chairs, beds, or other seating areas.

For added stability and comfort, the RS1 features padded lower leg supports that can be fully adjusted in width, angle, and height to suit each individual's needs. The ergonomically angled multi-grip handlebar is designed to support the natural "nose over toe" standing motion, encouraging proper positioning and balance during transfers. With its thoughtful design and adaptable support, the RS1 Patient Turner is a reliable solution for safer, more efficient patient mobility.

- Height adjustable handle
- Height adjustable knee pads
- Slip resistant turning disc
- Easy to use and manoeuvre
- Length: 472mm
- Width: 475mm
- Height: 914mm
- Castor Size: 75mm
- Max weight: 150kg
- BS EN ISO 10993-10



Max  
Weight  
150kg

90330



### RS2 Patient Transfer Assist

- Low step on height base for ease with clients who have minimal leg movement
- Anti slip base ensures the client is safe and secure in transfer
- Low friction castor for ease of manoeuvrability
- One foot touch brake keeps the base stable when embarking and disembarking from the raiser
- Versatile handle design ensures carers of all shapes and sizes can easily manoeuvre it
- Adjustable knee pad gives the client support and comfort when in use
- Length: 635mm, Width: 607mm, Height: 1186mm
- Max weight: 160kg
- BS EN ISO 10993-10

90329



### Ladder Belt (Sliding Support Strap)

- For use on RS2 Patient Transfer Assist
- Length: 2000mm
- Width: 100mm
- Height: 20mm

90359



### Patient Handling Belt

- For use on RS2 Patient Transfer
- 3 sizes available: Small, Medium and Large
- Length (mm): S:1050 / M:1300 / L:1500
- Width (mm): S:975 / M:1150 / L:1350
- Height: 100mm

Small	90333S
Medium	90333M
Large	90333L

## RS4 Patient Transfer Assist

- Light Weight, Strong and Robust frame
- Easy flip down support pads
- Easy to open leg mechanism. Works in conjunction with most types of seating and chairs
- Enables users to participate in transfers, maintaining patient independence
- Ideal for those users who can raise from a seated position to a standing position
- Adjustable kneepad
- 4" rear locking castors
- Compact narrow frame. Manoeuvre in small spaces
- Length: 820mm, Width: 530-730mm, Height: 1060mm
- Seat height: 620-670mm
- Max weight: 185kg
- BS EN ISO 10993-10



Semi Assembled

90724

# AQUA LIBRA

## STAINLESS STEEL SHOWER CHAIR RANGE



CONTACT OUR SALES TEAM ON  
**01249 760581**  
sales@urathon.com

# Aqua Libra AP

ATTENDANT PROPELLED

# Aqua Libra SP

SELF PROPELLED



3 HANDLE OPTIONS AVAILABLE



	Straight Handle		Push Handle		Pram Handle	
	AP	SP	AP	SP	AP	SP
<b>17"</b>	90646-50001	90649-50001	90646-50010	90649-50010	90646-50023	90649-50023
<b>19"</b>	90647-50002	90650-50002	90647-50011	90650-50011	90647-50024	90650-50024
<b>21"</b>	90648-50003	90651-50003	90648-50012	90651-50012	90648-50025	90651-50025



## FEATURES AND BENEFITS

- Fully modular construction – ability to adapt the chair to suit any user
- Stainless steel, white powder coated frame – cannot rust and lends itself well to constant recycling if required
- Height adjustment frame as standard, designed to be able to fit over most standard toilets
- Large commode pan capacity
- 17", 19" and 21" frame width options
- Max user weight: 160Kg
- Available in both attendant and self propelled variants
- Wide range of seat variances – ideal for all toileting requirements
- Adjustable / easily removable arms can be supplied with various width options
- Easily removable backrest with 3 handle options available
- Easy to clean seat and backrest materials
- Fully removable swing away footrests with simple length adjustment and nylon footplates.



### AP Aluminium Shower Commode Chair

- Removable PU Back and PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 460mm
- Seat height: 510-535mm
- Height: 985-1010mm
- Length: 860mm
- Width: 569mm
- BS EN ISO 10993-10



90209



### SP Aluminium Shower Commode Chair Wheel Kit

Turn your AP (90209) into a Self-propelled Aluminium Shower Commode Chair.

- 2 x Braked Nylon Castors and 2 x 24" MAG Wheels with PU Tyres and Q R pins
- Removeable PU back and horseshoe seat with aperture plug
- Fold up PU arm pads
- Removeable swing away leg rests and nylon angle adjustable footplates
- Potty
- Seat width: 460mm
- Seat height: 510-560mm
- Height: 990-1010mm
- Length: 1085mm
- Width: 660mm
- BS EN ISO 10993-10

90351



### AP HD Aluminium Shower Commode Chair

- Removable PU Back  
PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 560mm
- Seat height: 510-535mm
- Height: 985-1010mm
- Length: 860mm
- Width: 669mm
- BS EN ISO 10993-10



90350



## AP XXL Aluminium Shower Commode Chair

- Removable PU Back, PU Horseshoe Seat with Aperture Plug
- 4 x Braked Nylon Castors
- Fold Up PU Arms Pads
- Removable Swing Away Leg Rests and Nylon Angle Adjustable Footplates
- Potty
- Seat width: 600mm
- Seat height: 510-535mm
- Length: 860mm
- Width: 715mm
- Height: 985-1010mm
- Max weight: 300kg
- BS EN ISO 10993-10



90513



## AP Aluminium Shower Commode Chair, Tilt in Space

- Lightweight aluminium frame
- Easy to manoeuvre and offers long-term rust-free durability
- Tilt in space chair – ideal for showering comfortably
- 4 x Large, braked, robust castors for durability and manoeuvrability
- Polyurethane seat, seat-infill, back and armrests are water-resistant, easy to clean and comfortable
- Full height adjustable to suit all heights of the users
- Simple swing away and height adjustable footrests
- Easy to clean seat and removable back
- Flip up arm rests – For easy side transfer
- Large capacity commode pan and lid
- BS EN ISO 10993-10



90581

## Hospital Portering Chair

- Full PU Seat and back – Durable, robust, easy to clean and reduces any cross-infection risk. Ideal in heavy usage situations
- Strong balanced footplate – Durable with the ability to take the full user weight without tipping
- The footplate is made from HDPE plastic – Robust and ideal for infection control. Competitors tend to use wood which is a cross infection risk
- Stainless steel frame – Ensure product longevity. Allows the product to be used outside without fear of rust. Easy to wipe surface reduces risk of cross infection. Competitors are made from mild steel that can rust when exposed to the inclement weather
- 300kg SWL – Most competitors are only 190 KGs, so this can be used with the heaviest of patients
- Moulded drop-down arms – Allows for quick and safe side transfers
- Accessory mounts – Allows IV poles etc to be used with the chair
- Xray/notes pocket as standard – Enables the users records to be transported safely with the user from department to department.
- BS EN ISO 10993-10

- Size: 600mm x 865mm x 1040mm
- Seat Width: 540mm
- Seat Depth: 410mm
- Seat Height: 540mm
- Backrest Height: 480mm



90337



### Mild Steel Glideabout Commode

- 4 x Castors
- Removable, Grey Aperture Seat with Padded Vanity Cover
- Removable Arm Rests, Backrest, Leg Rests
- Potty
- Grey Vinyl Upholstery
- Length: 922mm
- Width: 530mm
- Height: 962mm
- Seat height: 528mm
- Max weight: 165kg
- BS EN ISO 10993-10

90517



### Height Adjustable Chrome Plated Glideabout Commode

- Height Adjustable
- Height: 985 - 1085mm
- Seat height: 525 - 625mm
- Max weight: 165kg

90294



### KD Glideabout Commode

- Easy to Dissassemble
- Removable Aperture
- Padded Vinyl Seat
- Removable Back, Arms, Swing Angle Leg Rests
- Length: 902mm
- Width: 520mm
- Height: 908mm
- Seat height: 530mm
- Max weight: 190kg
- BS EN ISO 10993-10

90519



### Shower Stool

- Height Adjustable
- Powder Coated Steel Frame
- Easy to remove PP Seat
- Aluminium option
- Depth: 370-395mm
- Width: 495-540mm
- Seat width: 380mm
- Seat depth: 280mm
- Seat height: 480-630mm
- Width between arms: 485mm
- Max weight: 160kg
- BS EN ISO 10993-10

90115



### Greenwall KD Ali Shower Chair

- Removable padded armrests
- Blow-molded seat and back provide comfort and strength and is warm to the touch
- Drainage holes in bench
- Aluminium frame
- Leg height adjusts in 25mm increments
- Weight: 4.4kg
- Seat width: 400mm
- Seat depth : 330mm
- Width Between Arms: 490mm

Max Weight 136kg

90666



**Toilet Surround Height and Width Adjustable**

- Height & Width Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 584 - 734mm
- Height: 645 - 795mm
- Max weight: 190kg
- BS EN ISO 10993-10

90438



**Toilet Surround**

- Height Adjustable
- Floor fixing option available
- Length: 450 - 455mm
- Width: 594mm
- Height: 645 - 795mm
- Max weight: 180kg
- BS EN ISO 10993-10

90702



90040  
**6"**

90041  
**8"**

90042  
**12"**

**Bath Seat**


- Strong suction feet
- Available in 150mm, 250mm and 350mm height options
- Length: 290mm
- Width: 450mm
- Max weight: 190kg



**4 Part Bath Step**

- Length: 375mm
- Width: 480mm
- Height: 110mm
- Max weight: 250kg

90116



**Male Urinal Bottle with lid**

- Length: 126mm
- Width: 98mm
- Height: 225mm

90597



**2"**

**Raised Toilet Seat, No Lid**

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90635




**4"**

**Raised Toilet Seat, No Lid**

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90636



**6"**

**Raised Toilet Seat, No Lid**

- Ergonomic design, Strong, durable and easy to clean
- Easy adjustable fixings for secure fit on most toilets
- Quick to fit or remove
- Max weight: 190kg

90637



**Static Shower Chair Height Adjustable**

- L: 480-505mm
- W: 620 - 645mm
- H: 775-925mm
- Seat H: 440-590mm
- Max weight: 160kg
- BS EN ISO 10993-10

90300



**Toilet Frame, Height Adjustable**

- Injection Moulded Seat
- L: 443-471mm
- W: 607-612mm
- H: 622-771mm
- Seat H: 428-578mm
- Max weight: 190kg
- BS EN ISO 10993-10

90341



**Swivelling Bath Chair**

- Adjustable Width
- Width: 770 - 820mm
- Height: 530mm
- Width between arms: 600-700mm
- Max weight: 150kg
- BS EN ISO 10993-10

90071



**Toilet Frame, Clip on Seat & Chute**

- Height Adjustable
- Length: 465 - 495mm
- Width: 558mm
- Height: 690 - 835mm
- Seat height: 483-628mm
- Max weight: 190kg
- BS EN ISO 10993-10

Version with lid available (90698)

90120



**26" Bath Board 4 Slat**

- Adjustable securing brackets
- Length: 240mm
- Width: 600mm
- Height: 20mm

90210



**27" Bath Board 4 Slat**

- Adjustable securing brackets
- Length: 240mm
- Width: 690mm
- Height: 20mm

90211



**28" Bath Board 4 Slat**

- Adjustable securing brackets
- Length: 240mm
- Width: 710mm
- Height: 20mm

90212



**26" Bath Board 6 Slat**

- Adjustable securing brackets
- Length: 360mm
- Width: 660mm
- Height: 20mm
- Max weight: 190kg

90254



**27" Bath Board 6 Slat**

- Adjustable securing brackets
- Length: 360mm
- Width: 690mm
- Height: 20mm
- Max weight: 190kg

90255



**28" Bath Board 6 Slat**

- Adjustable securing brackets
- Length: 360mm
- Width: 710mm
- Height: 20mm
- Max weight: 190kg

90256

# Crutch Caddy

GIVING INDEPENDENCE

2  
FREE  
RUBBER  
CRUTCH  
TIPS

An innovative mobility aid to safely store a pair of crutches in an upright position ready to use.



Easy to use  
lock and release



Suits all  
standard crutches\*



Glow in the  
dark

Home

Workplace

Hospitals

Waiting Rooms

Functions

Keeps crutches  
tidy and accessible



**2  
FREE  
RUBBER  
CANE  
TIPS**

# Cane Caddy

GIVING INDEPENDENCE



An innovative mobility aid to safely store a cane in an upright position ready to use.



Available in two hole sizes\*



Lightweight and portable

The unique Crutch and Cane Caddy may be positioned wherever the safe storage of crutches or canes are needed. By your bedside, in the bathroom or next to a comfortable chair. Lightweight and portable, they provide easy access to crutches or canes, giving independence and freeing up carer's time.

**Keeps canes tidy and accessible**



### Double Bedstick

- Adjustable Width
- Powder Coated White
- Length: 335mm
- Width: 960 - 1555mm
- Height: 440mm
- Max weight: 140kg
- BS EN ISO 10993-10

90354



### Double Bedstick

- Adjustable Width
- Chrome Plated
- Length: 335mm
- Width: 960 - 1555mm
- Height: 440mm
- Max weight: 140kg

90355



### Bedstick Slatted Bed

- Length: 822mm
- Width: 622mm
- Height: 500mm
- Max weight: 190kg
- BS EN ISO 10993-10

90321



### Bed Grab Rail

- Length: 653mm
- Width: 435mm
- Height: 510mm
- Max weight: 165kg
- BS EN ISO 10993-10

90521



### Bed Bumpers with Mesh Divide

- Blue with black mesh
- Length: 2000mm
- Width: 710mm
- Height: 10mm

90522




Max Weight 160kg

### Stacking Commode Chair

- Height Adjustable and Fixed Arms
- HF Vinyl Welded Padded Seat and Back
- Brown Powder Coated Steel Frame
- Anti Rattle Ring and DIP Moulded Ferrule
- Commode Pan and Lid
- L: 525-545mm W: 551mm H: 785-930mm
- Seat Height: 452-602mm
- BS EN ISO 10993-10

90096



Max Weight 160kg

### Stacking Commode Chair with Removable Arms

- Height Adjustable and Removable Arms
- HF Welded Padded Seat and Back
- Brown Padding and Brown Powder Coated Steel
- Anti Rattle Ring and DIP Moulded Ferrule
- Commode Pan and Lid
- L: 525-545mm W: 551mm H: 785-930mm
- Seat Height: 452-602mm
- BS EN ISO 10993-10

90323



### Overchair Table

This curved overchair/bed table is both width and height adjustable allowing the table to be set up perfectly to the requirements of the user. This table is great for use with chairs, single beds and wheelchairs.

The kidney shaped top is designed to curve gently around the user reducing the need for them to reach forward, it also allows them to rest their arms comfortably on the table top whilst still having a large usable work area.

- Height Adjustable
- Width: 1010mm
- Depth: 480mm
- Height: 710-910mm
- Max Weight: 5kg

90615



### Overbed Table

- Height Adjustable
- NO Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Top
- Length: 697mm
- Width: 413mm
- Height: 630 - 930mm
- BS EN ISO 10993-10

90290



### Overbed Table with Castors

- Height Adjustable
- 4 x 2" Braked Castors
- Cream Powder Coated Steel Frame
- Beech Veneer Table Top
- L: 697mm. W: 413mm
- Height: 700 - 1000mm
- BS EN ISO 10993-10

90291



### Perching Stool with Arms and HF Back

- HF Vinyl Seat
- Steel Loop Arms & HF Back
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 830 - 975mm
- Seat height: 540 - 685mm
- Max weight: 160kg
- BS EN ISO 10993-10

90245



### Perching Stool with Arms and Back

- HF Vinyl Seat
- Steel Loop Arms & Steel Back
- Adjustable Height
- Length: 415 - 460mm
- Width: 560 - 620mm
- Height: 815 - 960mm
- Seat height: 540 - 685mm
- Max weight: 160kg
- BS EN ISO 10993-10

90169



### Bariatric Perching Stool with PU Seat, PU Back & Arms

- Adjustable Height
- Length: 465-515mm
- Width: 650-710mm
- Height: 790-935mm
- Seat height: 485-630mm
- Max weight: 290kg
- BS EN ISO 10993-10

Max Weight 290kg

90346



### Perching Stool with PU Seat, Steel Arms and Back

- Adjustable Height
- Length: 415-460mm
- Width: 560-620mm
- Height: 815-960mm
- Seat height: 535-680mm
- Seat width: 325mm
- Seat depth: 265mm
- Max weight: 160kg
- BS EN ISO 10993-10

90272



### Perching Stool with PU Seat, Steel Arms and Padded Back

- Adjustable Height
- Length: 415-460mm
- Width: 560-620mm
- Height: 855-1000mm
- Seat height: 535-680mm
- Seat width: 325mm
- Seat depth: 265mm
- Max weight: 160kg
- BS EN ISO 10993-10

90273



### KD Perching Stool

- Removable Arms, Legs and Back
- Steel Arms
- Height Adjustable
- Max weight: 160kg
- BS EN ISO 10993-10

90626



### 2 Tier KD Dinner Trolley

- Height Adjustable
- 2 x Removable Shelves
- Length: 501mm
- Width: 445mm
- Height: 845 - 945mm
- BS EN ISO 10993-10

90072



### Half Step

- Adjustable Feet
- Removable slip resistant mat
- Plastic - will not rust
- Length: 430mm, Width: 735mm, Height: 112mm
- Max weight: 190kg

90141



### Two Tier Trolley with Brakes

- Push to lock brakes
- Length: 530mm
- Width: 595mm
- Height: 770 - 910mm
- Castor size: 150mm
- Max weight: 125kg
- BS EN ISO 10993-10

90515



### Back Rest

- Folds flat for easy storage
- Sturdy frame with breathable mesh fabric
- Adjustable soft pillow
- 5 levels of angle adjustment
- Width: 580mm, Depth: 515mm, Height: 560mm
- Max weight: 127kg

90668

## Universal Chair Raiser

Practical and cost-effective alternative to traditional riser recliner chairs, designed to assist users in achieving a safe standing position from their existing armchair. Engineered for use across community equipment services and home care settings, the unit can be fitted to most domestic armchairs or high-back chairs without requiring replacement of the user's preferred furniture.

- Remote Control
- Safe Working Load: 200kg
- Maximum User Weight: 160kg
- Adjustment: 3 positional height settings
- Installation Time: Approx. 10–15 minutes
- Application: Domestic armchairs and high-back chairs



90670

## Brown KD High Back Chair

The Urathon KD Brown High Back Chair has been designed as an adjustable height daily use arm chair, suitable for users up to 165kg (160kg when used in conjunction with universal chair raiser).

The chair features moulded arm rests, a comfortable, profiled, padded back, giving excellent lumbar support, easy to clean, flame retardent vinyl upholstery with optional wings for added comfort and head support.

Supplied in a carton flat-pack eliminating storage issues this chair when combined with the Urathon Universal chair raiser is a great solution for emergency hospital discharges that require swift provision of a riser chair.

**Universal Chair Raiser can be supplied prefitted to our Brown KD High Back Chair**



90764



90372

# yuwell Anytime™

CONTINUOUS GLUCOSE MONITORING SYSTEM

## A world-first rechargeable CGM transmitter

With readings every three minutes and unlimited data storage, this cutting-edge device uses Bluetooth technology to sync seamlessly with your phone or tablet. The Anytime CT-3 CGM transmitter lasts an impressive two years, reducing both costs and environmental waste. Its state-of-the-art ergonomic design ensures user-friendly operation, while superior sterile packaging enhances safety. Best of all, you only need to replace the sensor every two weeks - keeping the transmitter for long-term use!



NO FINGER  
PRICKS



NO SCANNING



COMFORTABLE  
TO WEAR FOR  
14 DAYS




CUSTOMISABLE  
ALERTS & ALARMS



RECHARGEABLE  
TRANSMITTER  
LASTS 2 YEARS

## Manage your glucose Anytime, Anywhere!




MARD  
6.9%

Prediabetic and Type 2

Order Now

[yuwellanytime.co.uk](https://yuwellanytime.co.uk)



**Urathon®**

# POWERCHAIR RANGE

PAGE 6

## HERO

Carbon fibre powerchair

Max Weight 135kg



PAGE 8

## NIRA

Dual battery powerchair

Range 13 miles



PAGE 7

## ZEUS

Auto-folding powerchair



PAGE 9

## AURA

Ultra-lightweight powerchair

Only 12.3kg



**5 YEARS WARRANTY**

Terms and Conditions apply

**DISCOVER MOBILITY FREEDOM**



# ***Urathon***<sup>®</sup>

## **UK WAREHOUSING**

Our dedicated UK warehouse allows for quick dispatch to you. No overseas shipping delays. Receive your products quickly, on time and when you need them. That's our promise to you.

## **SUPPLY FROM STOCK**

We fulfil product orders quickly from stock within our self managed UK warehouse. This allows us to always provide fast and reliable turnaround times with next day delivery available.

## **OUTSTANDING CUSTOMER CARE**

Customer care is at the heart of our business. Your satisfaction is paramount and we always aim to go above and beyond.

Call **01249 760581**

Email **sales@urathon.com**

Visit **www.urathon.com**

***Urathon***<sup>®</sup>

**SPECIALIST SUPPLIERS TO THE MOBILITY  
& REHABILITATION INDUSTRY**