



Introducing Urathon

Urathon have been designing and brand manufacturing for many of the healthcare industry's leading suppliers for almost 30 years, including many of the UK's largest providers to Nursing homes, Hospitals, NHS community equipment stores and Mobility retailers.

We offer a complete and seamless supply chain solution for all your Healthcare and Mobility product requirements, either from our large fully stocked warehousing in Wiltshire, UK or by full mixed containers under your own brand directly from our wholly-owned manufacturing base in China.

CIRRUS Rollator

Orange 90716 Grey 90717 Purple 90718

Lightweight aluminium rollator in three metallic colours. Easy fold frame, flat-free tyres, push-to-lock brakes, breathable fabric seat, one-touch handle height adjustment and removable shopping bag.

Only 7kg
Max weight 136kg



CIRRUS+ Rollator

90743

Lightweight rollator with 9" uprated flat-free wheels, padded seat and backrest, one-touch handle height adjustment and large removable zip-shut bag.

Max weight 136kg

Only 7kg



ALTO Rollator

Purple 90750 Orange 90751 Grey 90752

Lightweight modern 3-wheel rollator with removable storage bag, ergonomic handles with unique one touch height adjustment and self-locking easy-fold mechanism.

Only 5kg




NIMBUS Rollator

90695

One of the lightest carbon fibre rollators available. Flat-free tyres, padded seat, removable backrest, adjustable handles and push-to-lock brakes.

Only 6kg
Max weight 136kg

L: 720
W: 600
H: 800-980



Satin Carbon Fibre Finish



STRATUS Rollator

Pink 90734 Blue 90735

Bariatric Rollator with Suspension


Heavy-duty rollator supporting up to 227kg. Extra-wide padded seat, flat-free tyres, rear suspension, lightweight foldable frame and two stylish metallic finishes.



Max weight 227kg

Only 10kg

Seat Width (in) 24



ZYRA Rollator

90758

Lightweight 4-wheel rollator with three-position adjustable seat height, padded seat and backrest, adjustable handles, flat-free wheels and removable shopping bag.

Max weight 135kg

Only 7kg



LEO Rollator

Medium 90652 Small 90653

Easy-fold rollator with breathable seat, adjustable backrest, large flat-free wheels, adjustable handles and removable zip-shut storage bag.

Only 8kg

Max weight 130kg

L: 720
W: 590



ARIES Rollator

90654

Easy-fold rollator with breathable seat, adjustable backrest, large flat-free wheels, adjustable handles and removable zip-shut storage bag.

Only 7kg
Max weight 130kg

Seat Height: 620
Width (folded): 280



TAURUS Rollator

Medium 90656 Small 90657

Premium outdoor rollator with luxury stitched leather seat and bag, large treaded flat-free wheels, adjustable handles and padded backrest.

Max weight 130kg

Only 8kg



Stitched Leather



GEMINI Rollator

90703

Innovative dual-fold aluminium rollator with chunky flat-free tyres, padded seat, adjustable backrest, height-adjustable handles and removable shopping bag.

Only 8kg
Max weight 136kg

Folded L: 400
W: 235
H: 980



Urathon
SPECIALIST SUPPLIERS TO THE MOBILITY
& REHABILITATION INDUSTRY



NAIDEX 2026
BEST PRODUCT
IN SHOW

PRODUCT
SHOWCASE
2026



HALO FLAT Power Scooter 90762

The Halo Flat is a lightweight carbon fibre folding scooter. Compact design with adjustable steering column and digital display, padded seat with flip-up armrests, flat-free PU tyres and removable lithium battery.

Folds easily for transport and storage.



THERON Power Scooter 90760

Lightweight 16kg folding scooter with removable aircraft-safe battery, padded seat, fold-up armrests, simple controls, USB-C charging and up to 9.4 miles range per battery.




DRAGON Powerchair 90768

Innovative Power Wheelchair combining smart engineering with everyday practicality.

Advanced front omni wheels improve turning radius and manoeuvrability.

Also features adaptive seat-height system for seamless transitions to different social situations.

X-Fold technology ensures easy, space-saving storage.



HERO Powerchair Ruby 90744 Onyx 90756

Aviation grade carbon fibre powerchair with luxury padded seating, removable aircraft-safe lithium battery, fold-up armrests and powerful twin motors for hill climbing.



ZEUS Powerchair 90740

Stylish auto-folding powerchair with lightweight aluminium frame, removable aircraftsafe lithium battery, front suspension, padded upholstery and retractable carry handle.



NIRA Powerchair 90759

Lightweight aluminium folding powerchair with padded upholstery, removable aircraft-safe lithium battery, anti-tip bars, swing-away leg rests and front suspension.

A spare battery (included) is stored away and can be swapped during longer outings.



VELORA Powerchair 90755

Super-light carbon fibre folding powerchair with Bluetooth remote control, removable aircraft-safe lithium battery, unique length-adjustable and removable leg rests, padded upholstery, and anti-tip bars.



HELIA Powerchair 90765

Modern, lightweight powerchair that offers up to 9 miles range. An additional battery* can be stored in the rear of the chair.

Spacious 43cm luxury padded seat, carbon fibre frame and front wheel suspension.

* Not included as standard.



ORVIA Wheelchair 90681

Lightweight aluminium crash tested transit wheelchair with foldable frame, flip-up armrests, swing-away adjustable leg rests, padded upholstery and flat-free tyres.



APOLLO Wheelchair AP 90737 SP 90747

Super-light folding wheelchair with padded breathable upholstery, flat-free tyres, lap belt, anti-tip bars and swing-away adjustable leg rests.

Available in attendant propelled (AP) or self-propelled (SP) models.



YO-YO Wheelchair AP 90741 SP 90746

Lightweight folding wheelchair with full 4-wheel suspension, padded upholstery, adjustable armrests, swing-away leg rests and flat-free tyres.

Available in attendant (AP) or self-propelled (SP) models.



ZEPHYR Wheelchair AP 90742 SP 90748

Crash-tested lightweight aluminium wheelchair with padded upholstery, adjustable armrests, flat-free tyres and additional seat cushion.

Available in attendant propelled (AP) or self-propelled (SP) models.

